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# MultiTASC® Dissection/Transection Device

(Model Numbers 4500-03, 4500-04, 4500-05, 4500-06, 4500-07, 4500-08, 4500-09, 4500-10)

Instructions for Use - English

**STERILE EO**

## Description

The MultiTASC® Dissection/Transection Device consists of two coaxial metal rings which are connected to a removable handle through concentric tubes.

The device is designed to dissect an atheromatous core from the vessel wall. The core is then transected and extracted from the vessel.

## Indication for Use

This device is intended for remote endarterectomy of peripheral blood vessels during vascular reconstruction.

## Contraindications

The MultiTASC Device is not recommended for the following:

1. pregnant women,
2. patients with allergy to any of the device materials,
3. patients with an aneurysm in the targeted vessel,
4. patients with severe calcification,
5. use in coronary, internal carotid arteries, or other vessels of the central circulatory system.

## Warnings

1. Do not use if the package or device is damaged.
2. Do not use if inner packaging is opened outside a sterile environment.
3. This device is to be used by a qualified physician.
4. Do not use if vessel wall layers are overly coherent and the atheromatous segment is difficult to separate from the vessel.

## Precautions

1. United States Federal and other law restricts this device to sale on or by the order of a physician.
2. Do not use device past the expiration date printed on the labeling.
3. This device is single use only. Do not reuse, reprocess, or re-sterilize (reference "Re-sterilization/ Repackaging" section of this document for further details).
4. Never attempt to use this device without fluoroscopy.
5. Do not use excessive force to strip/remove atheromatous segments from occluded vessel.
6. Ensure ALL dissected atheromatous material is removed.
7. This device is not designed to transect calcified material. Attempted transection of such material may damage the device. If the rings or shaft of the device become bent, do not attempt to repair the device or use in the patient. Withdraw from the patient and dispose of the device.
8. Do not attempt dissection after core transection without inspecting the rings for misalignment. This may be done under fluoroscopy.

## Potential Complications

- Incisional/puncture site inflammation, hematoma, seroma, or bleeding
- Superficial or deep wound infection
- Vessel injury, dissection, perforation or rupture


- Vessel spasm
- Emboli, thrombosis or acute occlusion
- Fever, drug reactions or allergic reactions to contrast medium
- Additional surgical intervention
- Limb loss due to distal embolization
- Renal failure
- Re-occlusion/restenosis
- Post operative hemorrhage
- Edema
- Pseudo aneurysm
- Fractured atheromatous
- Death

#### Additional Required Items

- High resolution fluoroscopic imaging capabilities
- MultiTASC® Sizing Tool
- Radiopaque ruler/tape (optional)
- Periscope™ Dissector (as required)
- Martin™ Dissector (as required)
- EndoHelix™ Retrieval Device (as required)

#### Specifications

| Catalog Number | Description                                      | Inner Diameter | Outer Diameter | Useable Length | Sizing Tool |
|----------------|--|----------------|----------------|----------------|-------------|
| 4500-03        | MultiTASC Dissection/Transection Device; Size 5  | 3.5 mm         | 5.5 mm         | 55 cm          | 5           |
| 4500-04        | MultiTASC Dissection/Transection Device; Size 6  | 4.5 mm         | 6.5 mm         | 55 cm          | 6           |
| 4500-05        | MultiTASC Dissection/Transection Device; Size 7  | 5.5 mm         | 7.5 mm         | 55 cm          | 7           |
| 4500-06        | MultiTASC Dissection/Transection Device; Size 8  | 6.5 mm         | 8.5 mm         | 55 cm          | 8           |
| 4500-07        | MultiTASC Dissection/Transection Device; Size 9  | 7.5 mm         | 9.5 mm         | 55 cm          | 9           |
| 4500-08        | MultiTASC Dissection/Transection Device; Size 10 | 8.5 mm         | 10.5 mm        | 55 cm          | 10          |
| 4500-09        | MultiTASC Dissection/Transection Device; Size 11 | 9.5 mm         | 11.5 mm        | 55 cm          | 11          |
| 4500-10        | MultiTASC Dissection/Transection Device; Size 12 | 10.5 mm        | 12.5 mm        | 55 cm          | 12          |

### To Open Package


1. Open box and pull pouch out.
2. Peel apart pouch tabs and present product tray to personnel in sterile environment.

### Pre-Use Inspection

Prior to introducing device into the patient, complete the following functionality checks:

1. Inspect packaging and all components for damage.
2. Inspect rings for misalignment. If rings are misaligned, do not use device. A slight gap between the rings is acceptable.

Acceptable:

Unacceptable:

### Instructions for Use

1. Perform a surgical cut down to the vessel that requires treatment.
2. Isolate the vessel and perform an arteriotomy long enough to expose the atheromatous core.
3. Dissect the atheromatous core from the vessel wall using the appropriate surgical instruments.
4. Use the MultiTASC Sizing Tool to confirm the ring size required.
5. Select the MultiTASC Dissection/Transection Device that corresponds to ring size determined in the previous step.

6. Complete steps indicated in the "Pre-Use Inspection" section of this document prior to introducing device into the patient.
7. Device is packaged with torque handle in place. Slide torque handle to desired location and rotate torque handle cap to lock.

8. Place rings over the atheromatous core. Under fluoroscopic guidance, begin dissection by moving the device distally into the vessel. Slowly advance the device distally until the endpoint of atheromatous core is reached.


**Note:** The torque handle may be repositioned as the device is advanced. To do so, loosen the torque handle cap, reposition, and tighten torque handle cap to lock in desired position.

**Note:** A heavily calcified atheromatous core may require the use of the Martin Dissector. In the event that use of the Martin Dissector is required, remove MultiTASC Device and refer to the Instructions for Use accompanying the Martin Dissector.


9. When the endpoint of the atheromatous core is reached, loosen the torque handle cap completely and slide torque handle until entirely disengaged from the torque handle cap. Remove the torque handle from the shaft. Slide torque handle cap to retainer and tighten to secure.

10. Remove hub cap from proximal end of the shaft by rotating the hub cap.


**Note:** Do not rotate the hub or shaft while removing hub cap.

11. Verify that the actuator handle is in "Position 1" and slide over the hub into the mating slots.


12. Maintain the position of the hub/shaft and rotate the actuator handle in direction of arrow indicated on the actuator handle until a positive stop is detected (approximately 90°), locking hub/shaft and actuator handle together.


**Note:** Do not rotate the hubs or shaft while attaching the actuator handle.

13. To transect atheromatous core from the endpoint, rotate the t-handle from "Position 1" and retract proximally into "Position 2."

**Note:** Do not rotate actuator handle.

14. With atheromatous core transected from endpoint, advance the t-handle distally and rotate into "Position 3" to create a backstop for the atheromatous core for extraction.

15. Slowly retract entire device proximally to extract atheromatous core from vessel.

**Note:** Forceps may be required to guide the atheromatous core from vessel as it is being extracted.

Continue extraction until device is fully removed from the vessel. In the event that core fragments remain within the vessel, an EndoHelix Retrieval Device may be used. In the event that use of the EndoHelix Retrieval Device is required, remove the MultiTASC Device and refer to the EndoHelix Retrieval Device Instructions for Use.

16. Perform an angiogram to ensure the atheromatous core is completely removed from vessel and to determine the proper endovascular technique required to treat the intimal flap at the endpoint. If necessary, the Periscope Dissector Device may be used to aid in guidewire crossing of the endpoint. In the event that use of a Periscope Dissector Device is required, remove the MultiTASC Device and refer to the Periscope Dissector Instructions for Use.
17. In the event that further dissection is desired, remove actuator handle by rotating it in the direction opposite

the arrow indicated on the handle until a positive stop is detected (approximately 90°). Inspect the rings for misalignment, this may be done under fluoroscopy. If rings are misaligned, do not use device. Disconnect the actuator handle from the hub. Attach the hub cap to the hub. Loosen torque handle cap from the retainer and attach torque handle to the shaft. Position the torque handle in desired position and tighten torque handle cap. Return to Step 8 to complete dissection.

### **Storage**

Store in a cool, dry place.

### **Re-sterilization/Repackaging**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics and structural integrity of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

### **Safe Handling and Disposal**

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

### **Cleaning:**

1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

### **Packaging:**

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803

**Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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©2020 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH