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INFORMATION FOR USE

PRODUCT DESCRIPTION

The Wovex® prosthesis is a vascular prosthesis woven with PET (Polyethylene Terephthalate) yarn and impregnated with bovine origin collagen. This impregnation gives it a mean water permeability of less than $5 \text{ ml}\cdot\text{cm}^{-2}\cdot\text{min}^{-1}$, which minimizes bleeding during implantation, and eliminates the preclotting step. The bovine collagen is then gradually resorbed by the patient. The black guide lines and the crimped construction facilitate implantation. The Wovex® prosthesis is sterilized by irradiation.

INDICATIONS FOR USE

The Wovex® prosthesis is indicated for replacement or bypass procedures in aneurysmal diseases of the thoracic aorta, the abdominal aorta and visceral arteries. In these localizations, Wovex® prosthesis is also indicated for arterial reconstruction in patients requiring systemic heparinization.

CONTRAINDICATIONS

1. The Wovex® prosthesis is contraindicated for use in coronary arteries.
2. The Wovex® prosthesis is contraindicated in patients with known or suspected hypersensitivity to bovine collagen or polyethylene terephthalate.

WARNINGS

1. **DO NOT PRECLOT.**
2. **This is a single use device. Reuse, resterilization, reprocessing and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient. Any remaining part of the prosthesis after the procedure must be destroyed.**
3. **Do not use if the integrity of the packaging or the outer pouch or the sterility of the product is compromised (the outer pouch is a moisture barrier).**
4. **After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.**
5. **Use of Wovex® prosthesis in pregnant/breastfeeding women and pediatric patients has not been studied.**

PRECAUTIONS

1. The Wovex® prosthesis should only be used by surgeons qualified in vascular surgery and in an appropriate operative environment.
2. The outer pouch is a moisture barrier. If the outer pouch has been opened, destroy the prosthesis if it has not been used during the intended procedure.
3. During the clamping step, use only atraumatic, shod-type vascular clamps to avoid any damage to the prosthesis wall. Avoid repeated or excessive clamping in the same place on the prosthesis.
4. It is recommended that only taper point needles be used to perform anastomoses to avoid cutting the prosthesis fibers.
5. To reduce the risk of fraying, care should be taken when cutting woven prostheses. It is possible to use a cautery.
6. If a cautery is used, it is recommended to use a fixed temperature cautery.
7. For abdominal use, the lymphatics should be carefully ligated and sealed to minimize lymph collection around the prosthesis.
8. For abdominal use, care should be taken to separate the prosthesis and anastomoses from adjacent bowel structures with an adequate layer of retroperitoneal tissue or omentum/peritoneum to reduce the risk of secondary fistula formation or bowel erosion.
9. It is recommended that the prosthesis be handled with caution and that any excessive stretching is avoided, so as not to damage the collagen impregnation on the prosthesis.
10. In the event of leaks or any oozing of blood at the anastomoses or along the prosthesis, refer to the HEMOSTASIS NOTES section, in this information for use.

11. A date and a revision number for this information for use are included, for users' information, on the first page, below the LeMaitre Cardial address. If 3 years have elapsed between this date and the use of the product, the user should contact LeMaitre Cardial to find out if additional information is available.

ADVERSE REACTIONS

Adverse reactions that may occur with the use of this product or any vascular prosthesis include, but are not limited to, perioperative hemorrhage, thrombosis, thromboembolic complications, infection, anastomotic aneurysms, seroma formation, and aortoenteric fistula formation.

OPENING THE PACKAGE

The Wovex® prosthesis is packaged in double trays sealed with lids and the whole is inserted into a pouch. The prosthesis is sterile unless the inner package is opened or damaged. Open the package, observing the usual asepsis rules. Open the pouch, take the outer tray in one hand and peel back the lid. The inner tray can be handled using sterile gloves or instruments.

HEMOSTASIS NOTES

As with any vascular prosthesis, occasional difficulties with hemostasis may occur. In the event that hemostasis is not easily obtained, the surgeon may wish to consider the following :

- Investigate for systemic coagulopathy and treat appropriately.
- Compression, as necessary.
- Additional sutures and/or pledgets at the anastomosis, as necessary.
- Reversal of heparin with protamine sulfate, as necessary.
- Utilization of topical coagulation therapy, e.g., thrombin, as necessary.

WARRANTY

LeMaitre Cardial warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of the first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in LeMaitre Cardial sole discretion or refunding the net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL THE OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL LEMAITRE CARDIAL BE LIABLE TO YOU FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

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LeMaitre Cardial GmbH
LeMaitre Cardial SAS