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INFORMATION FOR USE

PRODUCT DESCRIPTION

The Dialine® II prosthesis is a vascular prosthesis knitted with PET (Polyethylene Terephthalate) yarn and impregnated with bovine origin collagen. This impregnation gives it a mean water permeability of less than 10 ml·cm⁻²·min⁻¹, which reduces bleeding to a minimum during implantation, and eliminates the preclotting step. The bovine collagen is then gradually resorbed by the patient.

The black guide lines and the crimped construction facilitate implantation.

The Dialine® II prosthesis is sterilized by irradiation.

INDICATIONS FOR USE

The Dialine® II prosthesis is indicated for replacement or bypass procedures in aneurysmal and/or occlusive diseases of the abdominal aorta, visceral arteries and peripheral arteries. The Dialine® II prosthesis is also indicated for arterial reconstruction in patients requiring systemic heparinization.

CONTRAINDICATIONS

1. The Dialine® II prosthesis is contraindicated for use in the thoracic aorta and coronary arteries.
2. The Dialine® II prosthesis is contraindicated in patients with known or suspected hypersensitivity to bovine collagen or polyethylene terephthalate.

WARNINGS

1. **DO NOT PRECLOT.**
2. **This is a single use device. Reuse, resterilization, reprocessing and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient. Any remaining part of the prosthesis after the procedure must be destroyed.**
3. **Do not use if the packaging presents an integrity defect which can compromise the sterility of the product or if the integrity of the outer pouch is compromised (the outer pouch is a moisture barrier).**
4. **After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.**
5. **Use of Dialine® II prosthesis in pregnant/breastfeeding women and pediatric patients has not been studied.**

PRECAUTIONS

1. The Dialine® II prosthesis should only be used by surgeons qualified in vascular surgery and in an appropriate operative environment.
2. The prosthesis diameter selection should take account of the stated diameter under pressure.
3. The outer pouch is a moisture barrier. If the outer pouch has been opened, destroy the prosthesis if it has not been used during the intended procedure.
4. During the clamping step, use only atraumatic, shod-type vascular clamps to avoid any damage to the prosthesis wall. Avoid repeated or excessive clamping in the same place on the prosthesis.
5. It is recommended that only taper point needles be used to perform anastomoses to avoid cutting the prosthesis fibers.
6. To minimize lymph collection around the prosthesis, the lymphatics should be carefully ligated and sealed.
7. In the abdomen, care should be taken to separate the prosthesis and anastomoses from adjacent bowel structures with an adequate layer of retroperitoneal tissue or omentum/peritoneum to reduce the risk of secondary fistula formation or bowel erosion.
8. It is recommended that the prosthesis be handled with caution and that any excessive stretching is avoided, so as not to damage the collagen coating on the prosthesis.

9. In the event of slight leaks or any oozing of blood at the anastomoses or along the prosthesis, refer to the HEMOSTASIS NOTES section, in this information for use.
10. A date and a revision number for this information for use are included, for users' information, on the first page, below the LeMaitre Cardial address. If 3 years have elapsed between this date and the use of the product, the user should contact LeMaitre Cardial to find out if additional information is available.

ADVERSE REACTIONS

Adverse reactions that may occur with the use of this product or any vascular prosthesis include, but are not limited to, perioperative hemorrhage, thrombosis, thromboembolic complications, infection, anastomotic aneurysms, seroma formation, and aortoenteric fistula formation.

OPENING THE PACKAGE

The Dialine® II prosthesis is packaged in double trays sealed with lids and the whole is inserted into a pouch. The prosthesis is sterile unless the inner package is opened or damaged. Open the package, observing the usual asepsis rules. Open the pouch, take the outer tray in one hand and peel back the lid. The inner tray can be handled using sterile gloves or instruments.

HEMOSTASIS NOTES

As with any vascular prosthesis, occasional difficulties with hemostasis may occur. In the event that hemostasis is not easily obtained, the surgeon may wish to consider the following:

- Investigate for systemic coagulopathy and treat appropriately.
- Compression, as necessary.
- Additional sutures and/or pledgets at the anastomosis, as necessary.
- Reversal of heparin with protamine sulfate, as necessary.
- Utilization of topical coagulation therapy, e.g., thrombin, as necessary.

WARRANTY

LeMaitre Cardial warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of the first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in LeMaitre Cardial's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL THE OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL LEMAITRE CARDIAL BE LIABLE TO YOU FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

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LeMaitre Cardial GmbH
LeMaitre Cardial SAS