Latis[®] Silicone Over-the-Wire Graft Cleaning Catheter

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Bruksanvisning - Svenska

Latis[®] siliconen Over-the-Wire Graft Cleaning Catheter

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Latis® Silicone Over-the-Wire Graft Cleaning Catheter

(Model Numbers A4GW6)

Instructions for Use - English

Rx only **STERILE** **EO** 

Ref	Length (cm)	Guidewire Compatibility	Liquid Volume (cc)	CO ₂ Volume (cc)	Inflated Diameter (mm)
A4GW6	60 cm	0.035	0.5	1.0	9 (±1)

Description

The Latis® Silicone Over-the-Wire Graft Cleaning Catheter consists of a catheter body with a Syntel® silicone balloon encased in a protective mesh on the distal end and two standard luer connectors on the proximal end. One lumen is used for balloon inflation and the other is used for infusion of liquids, or insertion over a guidewire.

Indications For Use

The Latis Silicone Over-the-Wire Graft Cleaning Catheter is indicated for the removal of thromboemboli from vascular grafts and for occlusion and infusion of fluids into a graft.

Contraindications

The Latis® Silicone Over-the-Wire Graft Cleaning Catheter is contraindicated for endarterectomy procedures and for use as a vessel dilator.

Precautions

- The Latis Over-the-Wire Graft Cleaning Catheter should only be used by a physician trained in endovascular procedures and familiar with the medical literature regarding techniques, complications, and hazards associated with the treatment of diseases of vascular grafts
- Prior to use, confirm that the catheter size is appropriate for the vessel to be accessed during the procedure
- Maximum recommended inflation volume should not be exceeded (See SPECIFICATION TABLE)
- Balloon degradation and rupture may result from deposits in the vessel, excessive hardening or exposure to adverse environmental conditions
- If difficulty or resistance is encountered during use, the cause should be determined and corrected before proceeding. Failure to do so may result in damage to the vessel or graft

Warnings

- Do not use if the product or sterile package is damaged
- For single use only, do not resterilize or reuse
- Do not autoclave
- Do not use in conjunction with a power injector
- Air should not be used to inflate the balloon
- The OTW Latis Graft Cleaning Catheter should not be used in native vessels

Instructions For Use

If using CO₂, corrections to inflation volumes may be necessary due to CO₂ diffusion through the balloon wall. For optimal control, use the smallest syringe capable of holding the inflation volumes.

1. Visually inspect the Latis Over-the-Wire Graft Cleaning Catheter package to confirm package integrity and completeness. If the package is open or damaged, or if any components are missing, do not use.
2. Prior to use, confirm compatibility with vascular access cannula.
3. Prior to use, the catheter should be purged of all air with the medium to be used for the procedure.
4. With the balloon deflated and the mesh cover in the collapsed position, insert the purged catheter into the graft, advance the catheter to the desired position. The catheter may be placed over a 0.035-inch guidewire if required.

5. Upon positioning the catheter tip just beyond the embolus or thrombus, inflate the balloon with sterile medium. Inflation should continue until the balloon can be felt to have engaged the graft wall.
6. Withdraw the catheter and occlusive material back through the graft access site. Ensure that the balloon is deflated if necessary to allow withdrawal through the arteriotomy or access sheath.
7. If additional passes are required during the procedure, repeat steps 3 to 5 after confirming the integrity of the balloon and catheter.
8. To infuse liquids, first remove any guidewire that may be placed in the through lumen/guidewire port, attach the infusion syringe to the luer lock hub.

Potential Complications

Potential complications associated with the use of the Latis Graft Cleaning Catheters include but are not limited to: systemic infection, local hematomas, vessel dissection, perforation and rupture, hemorrhage, vessel thrombosis, distal emboli of blood clots or arteriosclerotic plaque, air embolus, aneurysms, vessel spasms, arteriovenous fistula formations, and balloon rupture or tip separation with fragmentation and distal embolization.

Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

How Supplied

This device is supplied sterile and for single patient use. Pyrogen testing has been performed on this device. This device is supplied non-pyrogenic.

Storage Conditions

Handle with care. Product should be stored in a clean, cool, dry area away from chemical fumes.

Safe Handling and Disposal

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed of according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
 - a) Sodium hypochlorite solution (500-600 mg/l), or
 - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
 - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
 - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary

packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.

3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and Latis are registered trademarks of LeMaitre Vascular, Inc.
©2020 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH