**AnastoClip<sup>®</sup> Universal Clip Remover**

English — Instructions for Use

**AnastoClip<sup>®</sup> Universalclip-Extraktor**

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**Pince ôte agrafes universelles AnastoClip<sup>®</sup>**

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English — Instructions for Use

**STERILE EO Rx only**


**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

## Indications

The AnastoClip Universal Clip Remover is only intended to remove implantable clips (when necessary) during general surgery, vascular surgery, and neurosurgery.

The AnastoClip Universal Clip Remover is available in two (2) lengths.

Model Number	Usable Length	Width at rest	Application
A.) 4001-00, 4001-04, e4001-00, e4001-04	4" (10cm)	>24mm (.9")	Shallow, wide access
B.) 4013-00, 4013-01, e4013-00, e4013-01	6" (15cm)	>18mm (.7")	Deep, narrow access

*NOTE: Usable Length on the Remover refers to the distance between the tips and the midpoint of the serrated region on the handle.*

## Instructions For Use

To remove a misplaced clip, place the AnastoClip Remover underneath and through the apex of the already placed clip. Squeeze the handles of the AnastoClip Remover together. This process will open the clip and allow for the clip to be removed from the tissue.

## WARNING

The AnastoClip Remover is provided STERILE and is intended for use in a SINGLE procedure only.

## Complications

- Infection
- Pain/swelling
- Tissue damage
- Blood leakage
- Hemorrhage
- Intimal dissection
- Anastomosis rupture
- Dehiscence
- Necrosis
- Nerve injury
- Intimal hyperplasia
- Cerebrospinal fluid (CSF) leak
- Neuropathic pain
- Spinal infarction

**Resterilization/ Re-use**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the devices may be compromised due to reprocessing or re-sterilization since the devices were only designed and tested for single use. The shelf life of this device is based on single use only.

**Notices: Limited Product Warranty; Limitation Of Remedies**

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and AnastoClip are registered trademarks of LeMaitre Vascular, Inc.
© 2017 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH