HYDRO LeMaitre® Valvulotome

Instructions for Use - English

HYDRO LeMaitre® Valvulotom

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Valvulotome HYDRO LeMaitre®

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Valvulotomo HYDRO LeMaitre®

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Valvulotomo LeMaitre® HYDRO

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Valvulótomo HYDRO LeMaitre®

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STERILE **EO** **Rx only**

Description

The HYDRO LeMaitre® Valvulotome is a device that cuts venous valves during vascular procedures such as in-situ peripheral bypass, non-reversed translocated bypass, coronary artery bypass, and arterio-venous fistula creation. The Centering Hoops keep the head of the HYDRO LeMaitre Valvulotome centered in the vessel and prevents the valve-cutting blades from damaging the vessel wall. The size of the Centering Hoops and Cutting Blades adjusts to the internal diameter of the vein as the HYDRO LeMaitre Valvulotome is being drawn through the vessel. The LeMills Valvulotome may be provided as an accessory device to cut the last valve set as needed. Kit includes LeMills Valvulotome.

Indication for Use

The HYDRO LeMaitre Valvulotome is intended for the treatment of vascular disorders and more particularly for excising or disrupting venous valves.

The LeMills Valvulotome is intended to cut venous valves.

Contraindications

1. endarterectomy procedures,
2. thrombolysis procedures,
3. vein stripping procedures,
4. embolectomy procedures, or
5. vessel dilation procedures.

Warnings

1. Do not use if the package or device is damaged.
2. Do not use if inner packaging is opened outside a sterile environment.
3. The HYDRO LeMaitre Valvulotome and LeMills Valvulotome are to be used by a qualified physician.
4. Do not insert the HYDRO LeMaitre Valvulotome into a vessel or extract from a vessel in the open position.
5. Do not open or close the HYDRO LeMaitre Valvulotome while in a coiled configuration.
6. The HYDRO LeMaitre Valvulotome must be flushed with saline or heparinized saline only.
7. Do not pass the HYDRO LeMaitre Valvulotome or LeMills Valvulotome through a vessel that has undergone synthetic grafting or contains implants.
8. Do not rotate the HYDRO LeMaitre Valvulotome in a vessel.
9. Do not use device for valvulotomy unless the target vein is fully distended by arterial blood flow (due to anastomosis of artery to vein) or by saline injection.

Precautions

1. United States Federal and other law restricts this device to sale on or by the order of a physician.
2. Do not use devices past the expiration date printed on the labeling.
3. These devices are single use only. Do not reuse, reprocess, or re-sterilize (reference "Re-sterilization/Repackaging" section of this document for further details).
4. Keep any dissected portion of the vein straight and avoid twisting the vein during valvulotomy.
5. Use caution when tracking the HYDRO LeMaitre Valvulotome in thrombophlebitic veins.
6. Do not advance the HYDRO LeMaitre Valvulotome with the blades in the open position.
7. Do not rotate the Green Handle while the HYDRO LeMaitre Valvulotome device is in open position.

Potential Complications

- Vessel wall perforation
- Intimal disruption
- Residual tributaries
- Entrapment of cutting blades in branches or anastomoses
- Inadequate valve disruption


- Vein Stricture
- Thrombus formation
- Post-procedure wound necrosis
- Vessel damage
- Vessel Occlusion/ Stenosis
- Hematoma
- Hemorrhage
- Infections
- Erythema
- Entrapment of cutting heads/blades

Additional Required Items


1. Radiopaque Tape
2. 3cc Syringe

Specifications

HYDRO LeMaitre Valvulotome:

LeMills Valvulotome:

Catalog Number	Useable Length	Closed Catheter Outer Diameter	Maximum Open Cutting Blade Diameter	Maximum Open Centering Hoop Diameter	LeMills Valvulotome
1009-00, e1009-00	98 cm	1.5 mm	6.0 mm	9.5 mm	Included
1009-00J	98 cm	1.5 mm	6.0 mm	9.5 mm	Not Included
1010-00, e1010-00	40 cm	1.5 mm	6.0 mm	9.5 mm	Not Included
1011-00, e1011-00	98 cm	1.5 mm	6.0 mm	9.5 mm	Included
1012-00, e1012-00	40 cm	1.5 mm	6.0 mm	9.5 mm	Not Included

To Open Package

1. Open box and remove sealed sterile unit.
2. Open sealed sterile unit and present contents to personnel in sterile environment.
3. Carefully remove lids, as applicable.

Caution: Remove product from tray with care to ensure centering hoops are not damaged.

4. Place in basin filled with sterile, heparinized saline.

Pre-Use Check

Inspect blades for damage and alignment prior to use.

Instructions for Use


1. Cut the first valve under direct vision using standard technique.
2. Perform the proximal anastomosis to distend the vein.
3. With the HYDRO LeMaitre Valvulotome in the open position, remove the Irrigation Port cap and flush with sterile heparinized saline. Replace the cap upon completion.
4. Sheath the Centering Hoops within the Protective Sheath by advancing the White Handle until the Cutting Blades are fully enclosed and a positive stop is detected.
5. Activate the hydrophilic coating of the device by:
 - a) placing the HYDRO LeMaitre Valvulotome into a basin with enough sterile saline to cover the device, or
 - b) by wiping the Protective Sheath with sterile saline.
6. Advance the device into the vein to the treatment area; 2 cm to 3 cm distal to the anastomosis.
7. Expose the Cutting Blades and Centering Hoops by maintaining the position of the Green Handle and retracting the White Handle until a positive stop is detected.

Note: *Failure to maintain the position of the Green Handle while exposing the Cutting Blades may cause damage to the anastomosis.*

8. Slowly retract the device in the open position to disrupt the valve sets.

Warning: *Do not rotate the HYDRO LeMaitre Valvulotome.*

9. Proceed with caution until the Green Safety Stripes appear. These stripes indicate the 2 cm end of the Protective Sheath. When the third (and final) Green Safety Stripe is visible, sheath the Cutting Blades and Centering Hoops into the Protective Sheath by maintaining the position of the Green Handle and advancing the White Handle.

Warning: *Failure to sheath the Cutting Blades and Centering Hoops prior to retracting the device from the vein may cause damage.*

Note: *If sheathing the Cutting Blades for vein extraction leaves the most distal valve set intact, use the LeMills Valvulotome to excise using standard technique.*

Note: *If the device is set aside for further use within the same procedure, it must be stored in the open position in a basin of heparinized saline. Flush the device as noted in step 3 prior to reintroduction.*

10. Pass the device through the vein an additional 1-2 times.
11. Steps 5 through 9 may be repeated as necessary.
12. Confirm the free flow of blood through the vein.

Storage

Store in a cool, dry place.

Resterilization/Re-use

These devices are single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR

RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular GmbH