**LeverEdge<sup>®</sup> Contrast Injector**

Instructions for Use - English

**LeverEdge<sup>®</sup>-Kontrastmittelinjektor**

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**STERILE EO**

## Description

35 cc syringe with ISO 594 female luer lock.

## Indication for Use

The LeverEdge® Contrast Injector is intended to be used to perfuse contrast medias into the vessels for angiographic procedures.

## Contraindication

The LeverEdge Contrast Injector is not designed, sold, or intended for use except as indicated.

## Warnings

Read instructions carefully. The attending clinician is solely responsible for the setup and use of this device. Observe aseptic technique in all tubing connections. Contents of unopened, undamaged package are sterile and nonpyrogenic. Inspect packaging prior to opening. If damaged, sterility may be compromised. Inspect for damage before use. If damaged, do not use. Do not over-tighten rigid connections. This device is designed for single use only. Do not resterilize or reuse. Dispose of this device according to hospital procedure for contaminated material.

## Cautions

Federal (USA) law restricts this device to sale by or on the order of a physician. Observe all additional warnings and cautions contained in the following directions for use. Proper sterile surgical technique is necessary to use this product. The product is only to be used by a physician qualified in angiographic procedures. When infusing contrast media, techniques to avoid air or thrombus embolization should be observed.

## Directions for Use

The following procedure is for information purposes only. Proper surgical procedures and techniques are necessarily the responsibility of the medical professional. Each clinician must determine the appropriate use of this device based upon his or her medical training and experience and the type of surgical procedure employed.

1. Load LeverEdge Contrast Injector with contrast media using syringe loading technique.
2. Purge the system of all air bubbles by pointing device upward and pressing handle to expel.
3. Attach to the delivery system.
4. Press down on handle while watching angiographic image.
5. If more contrast is needed, push harder until handle reaches bottom of travel.

## Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

## Safe Handling and Disposal

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.


This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

## Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

## Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be

labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.

3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803

**Notices: Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS. These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction. A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and LeverEdge are registered trademarks of LeMaitre Vascular, Inc.
©2019 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH