**XenoSure<sup>®</sup> Biologic Patch**  
Instructions for use - English

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Processed Bovine Pericardial Patch

Instructions for Use - English

**STERILE | A**

## Storage

The XenoSure® should be stored between 0°C - 40°C (32°F-104°F). Avoid locations where extreme temperature fluctuations may occur; for example, near steam or hot water pipes, air conditioning ducts, or in direct sunlight. REFRIGERATION IS NOT REQUIRED. FREEZING WILL SERIOUSLY DAMAGE THE XENOSURE AND RENDER IT UNFIT FOR USE.

## Device Description

The XenoSure consists of one piece of bovine pericardial tissue that has been selected for minimal tissue blemishes. The tissue is treated with a glutaraldehyde process which crosslinks the collagen fibers and minimizes antigenicity. The XenoSure is liquid chemical sterilized and packaged in a plastic jar containing sterile glutaraldehyde storage solution. The XenoSure is designed to repair the body's natural organs.

The XenoSure comes in the following sizes:

Model			Size (cm)			Model			Size (cm)		
0.6P8	e0.6P8	0.6x8				2P9	e2P9	2x9			
0.8P8	e0.8P8	0.8x8				4P4	e4P4	4x4			
1P6	e1P6	1x6				4P6	e4P6	4x6			
1P10	e1P10	1X10				6P8	e6P8	6x8			
1.5P10	e1.5P10	1.5X10				8P14	e8P14	8x14			
1P14	e1P14	1x14				10P16	e10P16	10x16			

## Indications For Use

The XenoSure is intended for use as a surgical patch material for cardiac and vascular reconstruction and repair, soft tissue deficiency repair and reinforcing the suture line during general surgical procedures.

## Contraindications

None Known.

## Potential Complications

- Restenosis
- Pseudoaneurysm formation
- Infection
- Thrombosis
- Calcification
- Fibrosis
- Vessel occlusion
- Patch rupture
- Dilatation
- Myocardial infarction
- Bleeding

- Cerebrospinal fluid leakage
- Stroke
- Death

## Warnings

The principal complications that have been reported for bovine pericardial tissue are fibrosis and infection. These complications are observed only in a small minority of patients after implantation of the bovine pericardial tissue.

## Precautions

All persons responsible for the handling and preparation of the XenoSure Biologic Patch must exercise utmost care to avoid damage to the XenoSure Biologic Patch tissue.

- **FOR SINGLE USE ONLY.** Do not reuse, reprocess, or resterilize. Reuse, reprocessing, and/or resterilization of the device and/or failure could cause patient injury, illness or death. Any unused pieces of XenoSure Biologic Patch must be discarded. Note product "Use By" date.
- **INSPECT** sealed sterile package before opening. If seal is broken, contents may not be sterile and may cause infection in the patient. **DO NOT USE.** Do not discard the product. Please contact your distributor for further instructions.
- **DO NOT** expose the device to temperatures below 0°C (32°F) or above 40°C (104°F). **FREEZING WILL SERIOUSLY DAMAGE THE XENOSURE BIOLOGIC PATCH AND RENDER IT UNFIT FOR USE. DO NOT STORE UNDER REFRIGERATION.**
- **RINSE** the device according to the "RINSE PROCEDURE" section of this booklet before using. The XenoSure Biologic Patch storage solution contains glutaraldehyde and may cause irritation of skin, eyes, nose and throat. **DO NOT BREATHE STORAGE SOLUTION VAPOR.** Avoid prolonged skin contact and immediately flush area with water. In case of contact with eyes, seek medical assistance immediately. The liquid chemical storage solution should be disposed according to hospital procedure.
- **DO NOT** handle the XenoSure Biologic Patch with traumatic instruments. This may damage the device.
- **DO NOT** use any XenoSure Biologic Patch that has been damaged. Device integrity may be compromised.
- **DO NOT** attempt to repair the XenoSure Biologic Patch. Should damage to the XenoSure Biologic Patch occur before implantation, replace the XenoSure Biologic Patch.
- **DO NOT** expose the XenoSure Biologic Patch to antibiotics and or antimycotics as their substances may alter the cross linking of the aldehyde cured tissue.
- **DO NOT** resterilize. Unused sections should be considered non-sterile and discarded.
- **DO NOT** expose the XenoSure Biologic Patch to steam, ethylene oxide, chemical or radiation (gamma/electron beam) sterilization. Damage may result!
- **DO NOT** use cutting suture needles or cutting point armed sutures. This may damage the device.
- **DO NOT** allow the patch tissue to dry out during handling.
- **DO NOT** use if the device is beyond the expiration date.

## Adverse Effects

The XenoSure is designed to repair the body's natural organs. Improper functioning of an implanted XenoSure produces symptoms identical to symptoms that arise from deficiencies in the natural organ. It is the responsibility of the implanting surgeon to inform the patient of the symptoms that indicate improper functioning of the XenoSure.

1. Complete heart block and right bundle branch block are known complications reported for procedures involving cardiac repair near the A-V conduction bundles.
2. Glutaraldehyde-fixed bovine pericardium may be subject to late attack by

the immune system with subsequent tissue deterioration. The benefits of use of the XenoSure must be weighed against the possible risk of late tissue deterioration.

3. Animal studies with bovine pericardium have reported calcification and histological signs of deterioration as an adverse reaction. Findings include phagocytosis with accompanying chronic inflammatory infiltrate at the interface between bovine pericardium and surrounding host tissue with focal degradation of implant collagen consistent with host vs. graft reaction.
4. Bovine pericardium used for pericardial closure has been associated with epicardial inflammatory reactions and adhesions of the patch to the heart. Pericardial adhesions may increase the difficulty of repeat sternotomy.

### **How Supplied**

One XenoSure is provided sterile and non-pyrogenic in a sealed container; DO NOT RESTERILIZE. The patch is stored in a sterile phosphate buffered saline solution containing 0.2% glutaraldehyde. Sterility is assured if the package is unopened and has an undamaged seal. Unused sections should be considered non-sterile and discarded.

### **Directions For Use**

Choose the required XenoSure model as appropriate for the type of procedure being performed. The XenoSure can be cut to a size appropriate for a given repair. XenoSure is for SINGLE USE ONLY.

### **Patch Preparation**

Surgical gloves must be thoroughly washed to remove all powder residues before handling the XenoSure.

Examine the information on the jar label to verify selection of the correct XenoSure size. Carefully inspect the entire container and tamper-evident seal for damage.

DO NOT USE THE XENOSURE IF THE JAR IS DAMAGED OR IF THE SEAL IS BROKEN. Do not discard the product. Please contact your distributor for further instructions.

### **Rinse Procedure**

Remove the tamper-evident outer plastic seal and unscrew the jar cap. The contents of the jar are sterile and must be handled aseptically to prevent contamination. The outside of the jar is not sterile and must not enter the sterile field.

From the jar, remove the XenoSure by grasping its corners with sterile, atraumatic forceps.

Once removed from the container, submerge the XenoSure in the sterile saline. Using the same forceps, gently agitate the XenoSure in the basin. Allow the XenoSure to remain in the rinse basin until required by the surgeon.

Model		Size (cm)	Rinse Procedure
0.6P8	e0.6P8	0.6x8	500ml for 2 minutes
0.8P8	e0.8P8	0.8x8	
1P6	e1P6	1x6	
1P10	e1P10	1x10	
1.5P10	e1.5P10	1.5x10	
1P14	e1P14	1x14	
2P9	e2P9	2x9	
4P4	e4P4	4x4	
4P6	e4P6	4x6	
6P8	e6P8	6x8	1000ml for 3 minutes
8P14	e8P14	8x14	
10P16	e10P16	10x16	

### Alternate Patch Size Rinse Procedure

For any patch size not listed in the table, use the following rinse instructions.

Patch length (cm) x Patch width (cm) = Patch area (cm<sup>2</sup>)

If the patch has an area less than or equal to 24cm<sup>2</sup>, the patch requires a rinse time of 2 minutes in 500ml of saline.

If the patch has an area greater than 24cm<sup>2</sup> and less than or equal to 160cm<sup>2</sup>, the patch requires a rinse time of 3 minutes in 1000ml of saline.

*Note: The patch area is based on one side of the patch. This calculation has been formulated intentionally. Please follow the instructions.*

### Implantation

Cut and/or trim the XenoSure to the desired shape. Any excess XenoSure material should be treated as biological waste and discarded according to hospital procedure. During implantation, irrigate the XenoSure tissue frequently with sterile physiologic saline to prevent drying.

### Surgical Technique

It is beyond the scope of this Instructions for Use booklet to instruct the surgeon in specific repair procedures. LeMaitre Vascular, Inc. assumes that any surgeon performing the above operations has received adequate training and is thoroughly familiar with the pertinent scientific literature.

### Clinical Experience

Bovine pericardial patches have been successfully used in intracardiac repair procedures for the past 20 years with excellent clinical results. Published data on such use may be found in medical literature.

At this time, the long term durability of fixed bovine pericardial tissue is unknown. The clinical advantages in using this material along with its known characteristics must be weighed against its current undetermined long term durability.

### Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERMS INCLUDE

LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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## Symbol Legend

		Distributed By:			<b>Rx only</b>	<b>Open Here</b>	
English	Symbol Legend	Distributed By	Temperature limitations, store between 0°C and 40°C	Do Not Use if the Product Sterilization Barrier or its Packaging is Compromised	Caution: U.S. Federal and other law restricts this device to sale by or on the order of a physician	Open Here	Consult instructions for use: <a href="https://eifu.lemaitre.com">https://eifu.lemaitre.com</a>
Français Canadien	Légende des symboles	Distribué par	Conserver entre 0°C et 40°C	Ne pas utiliser si la jarre est endommagée ou si le sceau est rompu	Conformément à la loi américaine, ce produit ne peut être vendu que par un médecin ou sur ordonnance d'un médecin	Ouvrir ici	Consulter le mode d'emploi: <a href="https://eifu.lemaitre.com">https://eifu.lemaitre.com</a>

# XenoSure<sup>®</sup> Biologic Patch

Distributed By:

LeMaitre Vascular, Inc.  
Customer Service:  
Tel: (781) 221-2266  
Fax: (781) 221-2223

LeMaitre Vascular GK  
1F Kubodera Twin Tower Bldg.  
2-9-4 Kudan-minami, Chiyoda-ku  
Tokyo 102-0074, Japan  
Tel: +81-(0)3-5215-5681

LeMaitre Vascular ULC  
5915 Airport Road, Suite 608  
Mississauga, Ontario  
Canada L4V 1T1  
Tel: +1 905 673-2266

EC REP

LeMaitre Vascular GmbH  
Otto-Volger-Str. 5a/b  
65843 Sulzbach/Ts., Germany  
Tel: +49-(0)6196-659230

LeMaitre Vascular, Inc.  
63 Second Avenue  
Burlington, MA 01803

LeMaitre and XenoSure are registered trademarks of LeMaitre Vascular, Inc.  
©2014 LeMaitre Vascular, Inc. All rights reserved.