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Omniflow® II Biosyntetická cievna protéza

Návod na použitie - Slovenčina

Omniflow® II Biosynthetic Vascular Prosthesis

Instructions For Use - English

STORED IN
50% ETHANOL

STERILE A

Description

Omniflow® II Vascular Prosthesis is constructed of polyester mesh within an ovine fibrocollagenous tissue matrix. The prosthesis is sterilised in a glutaraldehyde solution.

The prosthesis is supplied sterile and nonpyrogenic in a solution of 50% ethanol. The prosthesis remains sterile unless the primary package is opened or damaged.

The Omniflow II Straight Vascular Prosthesis is mounted on a glass mandrel contained in a glass tube. The mandrel design prevents the prosthesis slipping off the mandrel when it is removed from the glass tube. The diameter and minimum length of the prosthesis is specified on the label applied to the glass tube.

The Omniflow II Curved Vascular Prosthesis is contained in a sterile flexible inner bag within an outer bag. The diameter and minimum length of the prosthesis is specified on the label applied to the outer surface of the outer bag.

The Omniflow II Vascular Prosthesis is considered MR safe.

Indications for Use

The Omniflow II vascular prosthesis is indicated for the replacement, reconstruction, bypassing or patching of diseased vessels in patients suffering occlusive or aneurismal diseases, in trauma patients requiring vascular replacement or for patients requiring vascular access such as for haemodialysis.

The Omniflow II Straight Vascular Prosthesis is intended:

1. To bypass, replace or reconstruct diseased or injured blood vessels
2. To patch and repair peripheral vessels
3. For arteriovenous access when a straight configuration is required.

The Omniflow II Curved Vascular Prosthesis is intended for arteriovenous access when a looped configuration is required.

Warnings

1. Do NOT re-sterilise the Omniflow II prosthesis. It is supplied sterile and pyrogen free. Use the prosthesis immediately after opening the package and discard any unused portions.
2. Do NOT use the prosthesis if the primary package is damaged as sterility may be compromised.
3. Do NOT use the prosthesis if the glass mandrel is broken.
4. Do NOT use the prosthesis if it is not completely covered by the storage solution.
5. Do NOT attempt to reposition the prosthesis after removal of the tunnelling instrument.
6. Do NOT straighten the curved prosthesis during preparation or implantation, as this will cause disruption of the mesh tissue interface.
7. Do NOT use the straight prosthesis to fashion a looped arteriovenous access as this may cause kinking.
8. Do NOT pull, stretch, twist, squeeze or pinch the body of the prosthesis.
9. Do NOT use ablation techniques such as cutting balloons, laser, or radio frequency ablation with the Omniflow II prosthesis.
10. Do NOT attempt to dilate the prosthesis with balloon angioplasty or stenting procedures.
11. The Omniflow II prosthesis should only be implanted by trained surgeons.
12. The use of the Omniflow II prosthesis in the coronary artery has not been evaluated.

Technical Information/Precautions

1. Ensure the rinsing procedure has been performed to remove the storage solution prior to implanting the prosthesis. Failure to do so may cause occlusion. Keep the prosthesis moist with sterile physiological saline during the procedure.
2. The use of a hollow tunnelling instrument for the passage of the prosthesis is essential. Failure to do so may cause disruption to the biosynthetic material and lead to occlusion, dilatation or aneurysm formation. The inner diameter of the tunnel should be at least 3mm larger than the indicated inner diameter of the prosthesis.
3. Ensure that the prosthesis does not become twisted when passing through the tunnelling instrument as this may lead to occlusion.
4. Avoid cross clamping with metal instruments as this may damage the prosthesis and cause occlusion, dilatation or aneurysm formation. If clamping is necessary use only a-traumatic clamps and avoid repeated or excessive clamping in the same position on the prosthesis.
5. The prosthesis has minimal longitudinal elasticity. Ensure the prosthesis is cut to the correct length. If it is too short it may cause suture pullout with a risk of anastomotic aneurysm. If it is too long it may kink and cause occlusion.
6. Cut off the sections of the prosthesis which were clamped during rinsing. Ensure that the full wall thickness and a mesh eyelet are incorporated with each stitch when performing the anastomosis. Failure to do so may result in stitch pullout and anastomotic aneurysm formation.

- Do not implant Omniflow II into a site with an active infection unless the surgeon determines there is not a more suitable alternative for preventing amputation or death.
- When the prosthesis is used for arteriovenous access some redness and swelling may be present over the implant area for a few days following implant.
- Insufficient data is available on which to base any conclusions regarding the use of the Omniflow II vascular prosthesis for aortocoronary bypass procedures.
- Failure to heparinize the prosthesis (i.e. in the case of patients who cannot tolerate heparin) may result in a higher likelihood of thrombosis or occlusion post-implantation, the extent of which has not been established.
- Omniflow II cannot grow in diameter or length and thus should not be implanted into infants or children unless a plan for its replacement is established and that no other suitable alternative treatment option exists.

Thrombectomy

Prompt intervention when occlusion occurs has been found to restore prosthetic function. Therefore the patient should be advised to report to the clinician if symptoms return. When performing a thrombectomy choose an embolectomy catheter of the appropriate size and follow the catheter manufacturer's instructions for use.

To avoid damage to the wall and flow surface of the prosthesis:

- Do NOT exert undue pressure when inflating the balloon
- Do NOT exert undue force when withdrawing the catheter from the prosthesis.

Contraindications

The prosthesis should not be used in patients with a known hypersensitivity to ovine material or glutaraldehyde.

Potential Complications with the use of Vascular Prostheses

Complications may occur with the use or in conjunction with any vascular prosthesis and include but are not limited to: infection, thrombosis/occlusion, dilatation, leakage, aneurysm formation, pseudoaneurysm formation, suture pullout and adverse tissue responses. Ongoing patient monitoring is recommended.

Pre-implant preparation

Prepare a sterile basin with adequate sterile physiological saline to cover the prosthesis, a sterile 20 ml syringe, sterile heparin 5,000 I.U. per ml.

Removing the prostheses from the container

STRAIGHT VASCULAR PROSTHESIS

- Remove the seal and cap from the glass tube. Remove the silicone stopper with sterile forceps. The forceps and stopper are discarded. Fig. 1.
- Aseptically remove the prosthesis from the glass tube by grasping the end of the glass mandrel with sterile forceps and lifting it gently out of the glass tube. Fig. 2. The hook on the bottom of the mandrel will prevent the graft from slipping off the mandrel.
- Grasp the mandrel by the hook, turn it upside down and allow the prosthesis to slide into the basin of sterile saline. If the prosthesis does not slide off easily, gently push it down the mandrel. Do NOT pull the prosthesis.

CURVED VASCULAR PROSTHESIS

- Aseptically remove the sterile inner bag containing the prosthesis from the outer bag. Fig. 3. Only the inner bag is sterile.
- Cut the corner of the inner bag with sterile scissors and allow the storage solution to drain into a container. Fig. 4. Discard the container of storage solution.
- Remove the prosthesis from the inner bag and place it into the basin of sterile physiological saline.

Fig. 1

Fig. 2


Fig. 3


Fig. 4

Rinsing the prosthesis

STRAIGHT VASCULAR PROSTHESIS

1. Use a sterile syringe to flush the lumen of the prosthesis with 20 ml of sterile physiological saline. Fig. 5. Repeat at least twice.
2. Clamp one end of the prosthesis. Insert sterile heparin at a concentration of 5,000 I.U. per mL into the lumen of the prosthesis so that the entire inside surface of the graft is coated with the solution (50,000 I.U. in 10 mL is sufficient for prostheses 35 cm and longer and 25,000 I.U. in 5 mL is sufficient for those shorter).
3. Apply a clamp to the other end of the prosthesis.
4. Place the prosthesis back into the basin completely covered by the sterile saline until it is time for implantation. This should be a minimum of 10 minutes.

When it is time to implant the prosthesis, remove the clamps from the prosthesis and allow the heparin and saline to drain out.

CURVED VASCULAR PROSTHESIS

The curved prosthesis is rinsed in the same manner as described above for the straight prosthesis. The curved configuration is maintained throughout the procedure. Fig. 6. Do NOT straighten the prosthesis.

Fig. 5

Fig. 6


Fig. 7


Fig. 8

Implantation

Prepare the incisions in the normal fashion.

Use of a hollow tunnelling instrument.

STRAIGHT VASCULAR PROSTHESIS

Form a tunnel for the prosthesis using a metal or plastic hollow tunnelling instrument and pass the prosthesis through the instrument. Fig. 7.

CURVED VASCULAR PROSTHESIS

Marking the proposed graft layout on the skin with a sterile surgical marker may facilitate final positioning of the prosthesis. Maintain the curve of the graft by inserting 2 hollow tunnelling instruments and passing both straight sections of the graft through simultaneously from the distal to the proximal incision. Fig. 8. Ensure that the graft does not kink at the curve during placement.

Carefully remove the tunnelling instrument when the graft has been correctly placed.

Anastomoses

1. When fashioning the ends for the anastomoses it is essential to cut off those portions of the prosthesis that were clamped, as clamping can damage the tissue matrix. No special cutting technique is required given the biosynthetic material of the graft.
2. Monofilament polypropylene is the preferred suture material. During suturing evert the edges of the prosthesis and ensure the full wall thickness and a mesh eyelet are taken up with each stitch. Fig. 9.
3. Prior to completion of the second anastomosis flush the prosthesis with blood to ensure complete removal of the heparin.
4. An intra-operative angiogram may be performed at the time of the procedure to document function. Ensure that the artery rather than the prosthesis is used for the injection.

At the completion of the procedure the patient record labels supplied with the prosthesis should be affixed in the patient and

hospital records.


Fig. 9

Arteriovenous Access

1. It is preferable to allow the prosthesis to heal for a minimum of two weeks prior to commencement of puncturing.
2. Do NOT puncture the arteriovenous access repeatedly in the same site as this may lead to disruption of the wall of the prosthesis, haematoma formation or pseudoaneurysm formation. Rotation of the puncture sites along the length of the prosthesis is necessary.
3. Do NOT puncture the curve of the prosthesis.
4. Do NOT puncture within 3 cm of the anastomotic sites.
5. Strict adherence to aseptic technique is essential to minimise the risk of infection.
6. Moderate compression of the puncture site will achieve haemostasis.

Vascular Patch

Cut the tubular prosthesis longitudinally to fashion into a vascular patch. Ensure that the blood flow surface is not damaged. When suturing the patch into place, ensure that the full wall thickness and the mesh are incorporated into each stitch.

Patient Information

The biosynthetic nature of Omniflow II does not change advice given to patients regarding post-operative care of the prosthesis as compared to other typical vascular prostheses. As such, good practice should include the following:

For bypass patients, avoid any compression of the prosthesis either by overly tight wrapping or clothing or laying on the prosthesis. Patients should report any abnormal signs immediately to the treating surgeon to reduce the risk of complications.

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1. Check the prosthesis every day. Feel for the pulse (thrill) in the prosthesis.
2. Inspect daily for signs of infection e.g. swelling, redness and pus.
3. Inspect daily for bruising or pseudoaneurysm formation.
4. Keep the skin over the prosthesis clean to help prevent infection.
5. Do not wear tight clothing or jewelry or drape or carry heavy items over the prosthesis as these may restrict blood flow.
6. Do not sleep on the arm with the prosthesis as this may occlude the blood flow.
7. Do not allow blood pressure testing or use of intravenous drips or the taking of blood from the arm with the prosthesis.
8. Patients should report any abnormal signs immediately to the dialysis unit to reduce the risk of complications.

Storage

The prosthesis must be kept in its original packaging until used. It must be kept in a dust free and dry environment at room temperature.

Return to Bio Nova International Pty. Ltd.

PATIENT IMPLANT DATA FORM

The implant data form should be completed immediately after surgery by the surgeon and returned in the envelope provided.

EXPLANTED SPECIMENS OF PROSTHESIS

If the Omniflow II vascular prosthesis is removed from a patient for any reason it may be returned to Bio Nova International Pty. Ltd. for analysis. Rinse blood and contaminants from the prosthesis using sterile saline. Return the specimen in a sealed container of 10% formalin. Include the serial number, catalogue number and reason for removal. Contact your local Omniflow II representative who will make arrangements for the specimen to be returned.

Restereilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot

be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. The unused section of the Omniflow II prosthesis or the explanted prosthesis should be considered as biohazard and be disposed per hospital policy or local regulations.

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre is a registered trademark of LeMaitre Vascular, Inc., and Omniflow is a registered trademark of Bio Nova Holdings Pty Ltd.
©2017 LeMaitre Vascular, Inc. All rights reserved.
LeMaitre Vascular GmbH