# Cardial

# BAIRD

## **COLLE CHIRURGICALE CARDIAL**

**Surgical Glue**

**Colle chirurgicale**

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**Colla Chirurgica**

**Cola quirúrgica**

**Chirurgische lijm**

**Cola Cirúrgica**

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# **INFORMATION FOR USE**

## **PRODUCT DESCRIPTION**

The Surgical Glue is a ready to use bi-component glue comprised of an adhesive (gelatin and resorcinol) provided in a tube, and a polymerizing agent (formaldehyde and glutaraldehyde) provided in an amber glass vial. The product is supplied with a specially designed syringe tip for accurate application of the polymerizing agent.

The Surgical Glue is supplied sterile. The adhesive is sterilized by the chemical action of the resorcinol. The polymerizing agent is sterilized by the chemical action of the glutaraldehyde and formaldehyde. The packaging and accessories are sterilized with ethylene oxide.

## **INDICATIONS FOR USE**

The Surgical Glue is indicated for joining dissected vessel layers and reinforcing sutures in acute aortic dissections, and reinforcing sutures in cardiac and vascular surgery.

## **CONTRAINDICATIONS**

The Surgical Glue is contraindicated in patients with known or suspected hypersensitivity to porcine gelatin.

## **WARNINGS**

- 1. The Surgical Glue is less effective in the presence of moisture or blood. The area to be glued must be as dry as possible before applying the adhesive. Failure to dry the area to be glued may result in longer polymerization times. In extreme cases, it may also result in inadequate adhesion, embolization, or tissue necrosis.**
- 2. DO NOT allow direct contact between the polymerizing agent and living tissues. The polymerizing agent is toxic and will cause localized tissue necrosis. Extreme care should be taken when injecting the polymerizing agent into the adhesive. Should contact occur, the tissues should be carefully wiped to remove the polymerizing agent.**
- 3. Use of Surgical Glue in pregnant/breastfeeding women and pediatric patients has not been studied.**
- 4. Safety and effectiveness of Surgical Glue in cerebrovascular repairs and intra-luminal areas have not been established.**
- 5. This is a single use device. Reuse, resterilization, reprocessing and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.**
- 6. Do not use if the packaging presents an integrity defect which can compromise the sterility of the product.**
- 7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.**

## **PRECAUTIONS**

- 1. Only physicians qualified in the appropriate surgical techniques should use Surgical Glue.**
- 2. An issue or revision date and a revision letter for these instructions are included for the user's information. In the event that 36 months have elapsed between this date and product use, the user should contact Cardial to see if additional product information is available.**

## **ADVERSE REACTIONS**

Adverse reactions that may occur with the use of this product include, but are not limited to tissue necrosis and embolization.

## OPENING THE PACKAGE

The Surgical Glue is packaged in double trays sealed with Tyvek lids. Open the package, observing the usual asepsis rules. Take the outer tray in one hand and peel back the lid. The inner tray can be removed using sterile gloves or instruments.

## INSTRUCTIONS FOR USE

1. Place the adhesive tube in sterile hot water before the surgical procedure. The adhesive is easier to apply when its temperature is raised to about 40°C.
2. Apply the adhesive on the area to be glued in a relatively thick layer.
3. Inject the polymerizing agent in the adhesive with the syringe. Two or three droplets of the polymerizing agent are generally necessary to polymerize about one milliliter of the adhesive. **SEE WARNING #2.**
4. During the short period of time required for polymerization, pressure should be applied to the glued area to obtain a satisfactory result. After polymerization, the color of the glue turns whitish and its adhesive power increases.

## WARRANTY

Cardial warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Cardial's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL CARDIAL BE LIABLE TO YOU FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

## References:

1. J. Bachet et al.:Surgery for Acquired Heart Disease. Aortic dissection - Prevalence, cause, and results of late reoperations. The Journal of Thoracic and cardiovascular surgery. Volume 108, Number 2, August 1994.
2. J. Bachet, B. Goudot, G. Teodori, D. Brodaty, D. Dubois, Ph. De Lentdecker, and D. Guilmet: Surgery of type A acute aortic dissection with Gelatine- Resorcine-Formol biological glue: A twelve-year experience. The Journal of Cardiovascular Surgery. 31 (3):263-273, May-June 1990.
3. J. Bachet: Dissections aortiques et colles - Application de la colle G.R.F. à la chirurgie de la dissection aigue de l'aorte. TechnoCoeur, Number 3, Mars-Avril 1992

## TRADEMARKS

Bard is a registered trademark of C. R. Bard, Inc., or an affiliate.

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## **Références :**

1. J. Bachet et al. : Surgery for Acquired Heart Disease. Aortic dissection – Prevalence, cause, and results of late reoperations. The Journal of Thoracic and Cardiovascular Surgery. Volume 108, Number 2, August 1994.
2. J. Bachet, B. Goudot, G. Teodori, D. Brodaty, D. Dubois, Ph. De Lentdecker, and D. Guilmet : Surgery of type A acute aortic dissection with Gelatine- Resorcine-Formol biological glue : a twelve-year experience. The Journal of Cardiovascular Surgery. 31(3) : 263-273, May-June 1990.
3. J. Bachet : Dissections aortiques et colles – Application de la colle G.R.F. à la chirurgie de la dissection aiguë de l'aorte. TechnoCoeur, Numéro 3, Mars-Avril 1992.

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## Literatur :

1. J. Bachet et al.:Surgery for Acquired Heart Disease. Aortic dissection - Prevalence, cause, and results of late reoperations. The Journal of Thoracic and cardiovascular surgery. Band 108, Nummer 2, August 1994.
2. J. Bachet, B. Goudot, G. Teodori, D. Brodaty, D. Dubois, Ph. De Lentdecker, and D. Guilmet: Surgery of type A acute aortic dissection with Gelatine-Resorcine-Formol biological glue: A twelve-year experience. The Journal of Cardiovascular Surgery. 31 (3):263-273, Mai-Juni 1990.
3. J. Bachet: Dissections aortiques et colles - Application de la colle G.R.F. à la chirurgie de la dissection aigue de l'aorte. TechnoCoeur, Nummer 3, Mars-Avril 1992

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1. J. Bachet et al.: Surgery for Acquired Heart Disease. Aortic dissection - Prevalence, cause, and results of late reoperations. The Journal of Thoracic and cardiovascular surgery. Volume 108, Number 2, August 1994.
2. J. Bachet, B. Goudot, G. Teodori, D. Brodaty, D. Dubois, Ph. De Lentdecker, and D. Guilmet: Surgery of type A acute aortic dissection with Gelatine- Resorcine-Formol biological glue: A twelve-year experience. The Journal of Cardiovascular Surgery. 31 (3):263-273, May-June 1990.
3. J. Bachet: Dissections aortiques et colles - Application de la colle G.R.F. à la chirurgie de la dissection aigue de l'aorte. TechnoCoeur, Number 3, Mars-Avril 1992

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