**LifeSpan® ePTFE Vascular Graft**

Instructions for Use – English

**LifeSpan® ePTFE Gefäßprothese**

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**STERILE**  **Rx only** 

## Product Description

The LifeSpan® ePTFE reinforced Vascular Graft consists of an expanded polytetrafluoroethylene (ePTFE) base tube reinforced with a thin spirally wrapped ePTFE tape. The graft also contains an orientation line.

The LifeSpan reinforced expanded PTFE stepped vascular graft also consists of an expanded polytetrafluoroethylene (ePTFE) base tube reinforced with a thin spirally wrapped ePTFE tape and an orientation line. This graft is configured with a 4 mm inner diameter end and a 7 mm inner diameter end with a step near the 4 mm end.

The LifeSpan expanded PTFE externally supported vascular grafts also consist of an expanded polytetrafluoroethylene (ePTFE) base tube reinforced with a thin spirally wrapped ePTFE tape and an orientation line. These grafts also have an external PTFE monofilament support for increased crush and kink resistance.

## Indications

- The LifeSpan ePTFE vascular grafts are indicated for use as a vascular prosthesis only. The grafts are intended for bypass or reconstruction of diseased or occluded blood vessels, or for arteriovenous shunts for blood access.
- The physician must evaluate each alternative method of treatment, discuss the risks and benefits with each patient, and decide whether to use a prosthetic vascular graft based upon all available factors.
- Grafts with removable external monofilament support over the length of the graft are used in bypass or reconstruction of occluded or diseased vessels, where compression or kinking could jeopardize patency.
- Grafts with external monofilament support in the middle of the graft may be used for the creation of an arteriovenous shunt for blood access; however, the graft must not be cannulated in the area of the external monofilament support.
- Stepped grafts are used for the creation of arteriovenous shunts for blood access. Stepped configurations may reduce the risk of steal syndrome and high cardiac output.

## Contraindications

- The LifeSpan ePTFE vascular graft is neither indicated nor approved for use in coronary artery bypass procedures (See Reference 1). Do not use the LifeSpan ePTFE vascular graft as a patch.
- LifeSpan ePTFE externally supported vascular grafts with external support along the entire graft length are not recommended for vascular access applications.

## Warnings

### For Single Use Only

1. Radiation should never be used for sterilization as this may damage the graft.
2. Do not expose the graft to either an open flame or to a laser without first providing adequate ventilation, as temperatures in excess of 400°C (752°F) produce highly toxic decomposition products.
3. Do not preclot the graft or immerse the graft in any solution (i.e. alcohol, oil, or aqueous solutions) prior to implantation. Exposing the graft to these solutions will affect the hydrophobic properties of ePTFE, resulting in wetting and possible excessive seroma formation.
4. Do not force any solutions through the graft lumen when examining the anastomosis during surgery as this may affect its hydrophobic properties, resulting in possible seroma formation.
5. Because the graft is not longitudinally elastic, the graft must be long enough to allow for a full range of body movements and to eliminate stress on the anastomoses.
6. Cut the graft with a sharp surgical instrument to avoid damage to the thin spirally wrapped ePTFE tape. Failure to correctly cut the LifeSpan ePTFE vascular graft may result in aneurysmal dilatation or reduced suture retention strength.
7. Do not pull or peel the thin spirally wrapped ePTFE tape from any area of the graft.
8. Use only clean gloves and atraumatic clamps or guarded (i.e. rubber shod) clamps to prevent damage to the graft wall during implantation. Avoid repeated or excessive clamping on any section of the graft.
9. Protect graft from damage by heavy or sharp objects.
10. When using a graft with removable external PTFE monofilament support, slowly unwind the monofilament at a

right angle to the graft. Be careful to avoid removal of the thin spirally wrapped ePTFE tape. Rapid unwinding and/or removal of the monofilament support parallel to the axis of the graft may damage the product.

11. Externally supported grafts used for arteriovenous shunts for blood access must not be cannulated in the area of the external monofilament support.

## Surgical Precautions

### Cutting

- Before cutting the graft, gently pull the graft taut and determine the correct length. The patient's body weight and posture must be considered as well as the anatomical location of the graft. It is especially important that the graft be cut long enough for axillofemoral, axillobifemoral, femorofemoral bypass or any other extraanatomic procedures.
- After adequate control of the distal vessel is achieved, use a new knife blade to make a longitudinal incision in the artery at the site of anastomosis. Utilize sharp scissors to cut one end of the ePTFE graft at a bevel to obtain an opening equal to the length of the arteriotomy.
- When cutting the graft, a sharp surgical instrument must be used to prevent damage to the thin spirally wrapped ePTFE tape. If the thin spirally wrapped ePTFE tape becomes frayed after cutting the graft, trim that portion of graft with a sharp surgical instrument. The thin spirally wrapped ePTFE tape should never be pulled or peeled away from any portion of the graft. Removal of this outer layer could result in aneurysmal dilation or a reduction in suture retention strength.

### Tunneling

- Use a tunneler to create a tissue tunnel that closely approximates the graft diameter. An oversized tissue tunnel may result in delayed or insufficient perigraft tissue attachment and may promote perigraft seroma formation.
- Extreme care should be taken to prevent kinking of the graft when passing it through tunnels.
- When creating a tunnel for the LifeSpan ePTFE externally supported grafts, make the tunnel one to two millimeters larger than the internal diameter of the graft being implanted to allow free passage of the external monofilament support.

### Suture Technique

Disruption of the graft, host vessel, and/or suture lines may occur if appropriate suture techniques are not observed.

- Utilize a small diameter, tapered, noncutting needle to prevent bleeding from suture holes. A full-radius cutting needle may damage the graft and should not be used.
- During anastomosis, use appropriate suture placement and small equal bites and avoid undue tension on the suture line.
- Following completion of the distal anastomosis, the graft should be gently filled and flushed with heparinized saline solution until circulation is restored. Be careful not to exert too much pressure with solution as possible seroma formation may result.

## Special Considerations

### Extraanatomic Bypasses: Axillofemoral and Femorofemoral Grafts

The following techniques must be used when performing these procedures:

- Prepare and drape the patient to allow full extension of the arm, shoulder girdle, or legs when determining appropriate graft length (refer to Warning #6 and Surgical Precautions subsection: Cutting). Prep the arm circumferentially to facilitate its movement during the procedure to check tension of the graft and prevent kinking.
- Separate the fibers of the pectoralis major and divide the pectoralis minor tendon to avoid angulation of the graft. The first and second portions of the axillary artery must be dissected and controlled and appropriate area for anastomosis selected.
- When dissections are complete, the subcutaneous tunnel should be made by curving to the mid-axillary line to avoid angulation over the costal margin in the sitting position.
- Position the axillary anastomosis medially in the first third of the axillary artery in the subfascial plane under the pectoralis major muscle and the external oblique fascia to aid in preventing anastomotic disruption. The third portion of the artery should never be used for bypass because shoulder motion affects this segment.
- Anastomose the graft as perpendicular to the axillary artery as possible to minimize stress. The anastomotic angle should be as small as possible and should not exceed 25° relative to the cut edge of the graft.

- Avoid protracted hyperabduction of the arm. Prolonged hyperabduction may lead to brachial plexus injury.
- To promote satisfactory healing, the patient should avoid extreme or abrupt movements of the arm, shoulder, or legs during a post-operative period of 1.5 to 2 months. Specifically, the patient should not reach out in front, raise arms above shoulder level, throw, pull, stride, or twist.

### **Vascular (Blood) Access**

- Patient monitoring is essential when the graft is used for vascular access to prevent excessive damage from complications such as infiltration, infection, thrombosis, embolism etc.
- To minimize subcutaneous bleeding during implantation, create a subcutaneous tunnel that approximates the size of the diameter of the graft.
- Rotate cannulation sites to prevent complications such as a disruption of the graft material and formation of a perigraft hematoma or pseudoaneurysm.

### **Post Surgery Invasive Procedures**

- LifeSpan ePTFE grafts are not radially elastic. Therefore, when performing an embolectomy using balloon angioplasty catheters within the lumen of the graft, the inflated balloon size must match the inner diameter of the graft. Over-inflation of the balloon or use of an inappropriate sized catheter may damage the graft and/or the balloon.
- If thrombectomy is necessary during the early post-operative period, a longitudinal incision with stay sutures is recommended.

## **Possible Complications with Vascular Protheses**

Failure to follow proper techniques discussed under Warnings, Surgical Precautions, and Special Considerations may result in:

- Mechanical disruption or tearing of the graft, suture line, or host vessel.
- Extreme blood loss
- Loss of limb function, limb, or death.

In addition, life threatening complications which may occur in conjunction with the use of any vascular prosthesis include, but are not limited to: Excessive suture hole bleeding; thrombosis; thromboembolic complications; infection; ultrafiltration or perigraft seroma; swelling of limbs; pseudoaneurysms; perigraft hematomas; skin erosion; steal syndrome; preoperative hemorrhage; aortoenteric fistula.

Complications due to individual patient reaction to an implanted device, or to physical or chemical changes in the components, may necessitate reoperation and replacement (sometimes within hours or days) of the prosthetic device. Careful and continuous medical follow-up is advised so that prosthesis-related complications can be diagnosed and properly managed to minimize danger to the patient.

## **Sterility**

The LifeSpan ePTFE vascular graft is supplied sterile and nonpyrogenic if the package has not been previously opened and is undamaged. Provided that the integrity of the package is not compromised, the package will serve as an effective barrier for a minimum of 5 years from the sterilization date. The graft is supplied in double plastic trays to facilitate handling and transfer to the sterile field during surgery. After opening the outer tray, the inner tray may be placed directly into the sterile field.

To open, hold the base of the outer tray, peel back the lid so that the inner tray can be removed by grasping the sealed lip. Beginning at one corner, peel back the inner tray lid and gently remove the graft. Use clean gloves or atraumatic instruments when handling the graft.

## **Resterilization**

When resterilization is performed, each institution should establish sterilization procedures that include biological indicators to establish the efficacy of their procedure. Sterilizers should be validated on a regular basis. The LifeSpan ePTFE vascular grafts can be resterilized ten times by steam as long as contamination with blood or any foreign material has not occurred. To resterilize, gently remove the graft from the plastic tray with clean atraumatic instruments or gloved hands and place in a hospital sterilization pouch. Do not use tray to perform autoclave resterilization procedures. Do not sterilize with any type of radiation sterilization equipment.

Normal autoclave cycle (gravity displacement steam sterilizer):

- 250°F (121°C) for 30 - 40 minutes or
- 270°F (132°C) for 15 - 20 minutes.

Flash autoclave cycle (prevacuum steam sterilizer):

270°F (132°C) for 4 minutes.

LeMaitre Vascular, Inc. does not guarantee the sterility of grafts after resterilization.

LeMaitre Vascular, Inc. has not qualified the use of ethylene oxide sterilization with this graft, and therefore; cannot provide recommendations for resterilization of the graft using ethylene oxide gas.

## Storage

To minimize contamination and to provide maximum protection, the prosthesis in the outer wrap should be stored in the outer cardboard box in a dry, contamination-free area until needed. Stock rotation at regular intervals is recommended to ensure usage of the graft prior to the “use by” date stamped on the label.

## Case History

### *Recovered Implants*

LeMaitre Vascular, Inc. is extremely interested in obtaining recovered clinical specimens of LifeSpan ePTFE vascular grafts for analysis. A written report summarizing our findings will be provided upon completion of our evaluation. Please contact your local company representative for return of recovered grafts. The grafts should be placed in a solution of 10% formalin or 2% glutaraldehyde. Refrigeration is not necessary under these circumstances.

## Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERMS INCLUDE LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular’s sole option) following the purchaser’s return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user’s information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

## References

1. Ascer, Enrico, M.D., Sushil K. Gupta M.D., Frank J. Veith, M.D., and Kurt R. Wengerter, M.D. Expanded Polytetrafluoroethylene Vascular Grafts, in *Vascular Grafts*, W. B. Saunders, 1989, pp. 460-473.
2. Brener, Bruce J. and Donald K. Brief. *Extraanatomic Bypasses in Vascular Surgery Principles and Practice*, McGraw-Hill Book Company, 1987, pp. 414-424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. *Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts*, in *Vascular Graft Update: Safety and Performance*, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82-94.
4. Johnston, K. Wayne, M.D., F.R.C.S. *Prevention and Management of Common Complications of Vascular Surgery—An Overview*, in *Vascular Grafts*, W. B. Saunders, 1989, pp. 487-491.

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### **Literaturhinweise**

1. Ascer, Enrico, M.D., Sushil K. Gupta M.D., Frank J. Veith, M.D., and Kurt R. Wengerter, M.D. Expanded Polytetrafluoroethylene Vascular Grafts, in *Vascular Grafts*, W. B. Saunders, 1989, pp. 460-473.
2. Brener, Bruce J. and Donald K. Brief. Extraanatomic Bypasses in *Vascular Surgery Principles and Practice*, McGraw-Hill Book Company, 1987, pp. 414-424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts, in *Vascular Graft Update: Safety and Performance*, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82-94.
4. Johnston, K. Wayne, M.D., F.R.C.S. Prevention and Management of Common Complications of Vascular Surgery—An Overview, in *Vascular Grafts*, W.B. Saunders, 1989, pp. 487-491.

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## Références

1. Ascer, Enrico, M.D., Sushil K. Gupta M.D., Frank J. Veith, M.D., and Kurt R. Wengerter, M.D. Expanded Polytetrafluoroethylene Vascular Grafts, in *Vascular Grafts*, W. B. Saunders, 1989, pp. 460-473.
2. Brener, Bruce J. and Donald K. Brief. *Extraanatomic Bypasses in Vascular Surgery Principles and Practice*, McGraw-Hill Book Company, 1987, pp. 414-424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. *Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts*, in *Vascular Graft Update: Safety and Performance*, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82-94.
4. Johnston, K. Wayne, M.D., F.R.C.S. *Prevention and Management of Common Complications of Vascular Surgery—An Overview*, in *Vascular Grafts*, W.B. Saunders, 1989, pp. 487-491.

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2. Brener, Bruce J. and Donald K. Brief. *Extraanatomic Bypasses in Vascular Surgery Principles and Practice*, McGraw-Hill Book Company, 1987, pp. 414-424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. *Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts*, in *Vascular Graft Update: Safety and Performance*, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82-94.
4. Johnston, K. Wayne, M.D., F.R.C.S. *Prevention and Management of Common Complications of Vascular Surgery—An Overview*, in *Vascular Grafts*, W.B. Saunders, 1989, pp. 487-491.

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## Referencer

1. Ascer, Enrico, M.D., Sushil K. Gupta M.D., Frank J. Veith, M.D., and Kurt R. Wengerter, M.D. Expanded Polytetrafluoroethylene Vascular Grafts, in Vascular Grafts, W. B. Saunders, 1989, pp. 460-473.
2. Brener, Bruce J. and Donald K. Brief. Extraanatomic Bypasses in Vascular Surgery Principles and Practice, McGraw-Hill Book Company, 1987, pp. 414-424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts, in Vascular Graft Update: Safety and Performance, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82-94.
4. Johnston, K. Wayne, M.D., F.R.C.S. Prevention and Management of Common Complications of Vascular Surgery—An Overview, in Vascular Grafts, W.B. Saunders, 1989, pp. 487-491.

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## Referenser

1. Ascer, Enrico, M.D., Sushil K. Gupta M.D., Frank J. Veith, M.D., and Kurt R. Wengerter, M.D. Expanded Polytetrafluoroethylene Vascular Grafts, in Vascular Grafts, W. B. Saunders, 1989, pp. 460–473.
2. Brener, Bruce J. and Donald K. Brief. Extraanatomic Bypasses in Vascular Surgery Principles and Practice, McGraw-Hill Book Company, 1987, pp. 414–424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts, in Vascular Graft Update: Safety and Performance, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82–94.
4. Johnston, K. Wayne, M.D., F.R.C.S. Prevention and Management of Common Complications of Vascular Surgery—An Overview, in Vascular Grafts, W.B. Saunders, 1989, pp. 487–491.

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## References

1. Ascer, Enrico, M.D., Sushil K. Gupta M.D., Frank J. Veith, M.D., and Kurt R. Wengerter, M.D. Expanded Polytetrafluoroethylene Vascular Grafts, in Vascular Grafts, W. B. Saunders, 1989, pp. 460-473.
2. Brener, Bruce J. and Donald K. Brief. Extraanatomic Bypasses in Vascular Surgery Principles and Practice, McGraw-Hill Book Company, 1987, pp. 414-424.
3. Colone, William M., Michael E. McClurken, Ph.D., and James M. McHaney. Physical Properties and Test Methods for Expanded Polytetrafluoroethylene (PTFE) Grafts, in Vascular Graft Update: Safety and Performance, sponsored by ASTM Committee F-4 on Medical and Surgical Materials and Devices and Subcommittee F04.04 on Cardiovascular Implants and Materials; Cosponsored by Food and Drug Administration, International Center for Artificial Organs and Transplantation, and Department of Health and Human Services, published by American Society for Testing and Materials, 1986, pp. 82-94.
4. Johnston, K. Wayne, M.D., F.R.C.S. Prevention and Management of Common Complications of Vascular Surgery—An Overview, in Vascular Grafts, W. B. Saunders, 1989, pp. 487-491.

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LeMaitre Vascular, Inc.  
63 Second Avenue  
Burlington, MA 01803

LeMaitre and LifeSpan are registered trademarks of LeMaitre Vascular, Inc.

©2011 LeMaitre Vascular, Inc. All rights reserved.

R2270-01 Rev. A 09/11

0088

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LeMaitre Vascular GmbH