Cardial

BAIRD

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Valvulotome

Dévalvuleur

Valvulotome

Valvulotomo

Valvulotomo

Valvulotoom

Valvulótomo

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INFORMATION FOR USE

PRODUCT DESCRIPTION

The Chevalier valvulotome consist of a 70 cm stainless steel shaft with two tips :

- The cutting tip is shaped like an inverted tulip with four teeth at its base. This tip is designed to cut the cusps of the valves.
- The other tip is atraumatic and used to test the patency and the course of the vein before the valvulotomy.

The Chevalier valvulotome is available in three tip diameters: 2.5 mm, 3.0 mm and 3.5 mm.

It is supplied sterile.

INDICATIONS FOR USE

The Chevalier valvulotome is designed to perform the valvulotomy of veins which are to be used as grafts in the in situ or non reversed bypass procedure.

CONTRAINDICATIONS

No contraindications are known for the Chevalier valvulotome.

WARNINGS

1. **The Chevalier valvulotome should not be passed through the proximal anastomosis. If the proximal valves are close to the long saphenous vein arch disinsertion, they can be cut with scissors under visual control.**
2. **The Chevalier valvulotome should be handled with care during the procedure to avoid any vein wall disruption or puncture.**
3. **This is a single use device. Reuse, reesterilization, reprocessing and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.**
4. **Do not use if the packaging presents an integrity defect which can compromise the sterility of the product.**
5. **After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.**

PRECAUTIONS

1. Only physicians qualified in the appropriate surgical techniques should use the Chevalier valvulotome.
2. An issue or revision date and a revision letter for these instructions are included for the user's information. In the event that 36 months have elapsed between this date and product use, the user should contact Cardial to see if additional product information is available.

ADVERSE REACTIONS

No adverse reactions have been reported for this device.

OPENING THE PACKAGE

The Chevalier valvulotome is packaged in double sealed pouches. Open the package observing the standard aseptic technique. Take the outer pouch in one hand and peel back. The inner pouch can be removed using sterile gloves or instruments.

INSTRUCTIONS FOR USE

1. The Chevalier valvulotome has been designed to perform valvulotomy of veins which are primed.
2. If necessary, the patency and the course of the vein should be tested with the atraumatic tip, before the valvulotomy.
3. After the proximal anastomosis has been performed, the correct size cutting tip is inserted distally and advanced proximally. It should not be passed through the proximal anastomosis.
4. The valvulotomy is performed by withdrawing the valvulotome back through the filled valve where a small resistance may be felt. The instrument is then passed back and forth through the valve a few more times, usually four or five, rotating a quarter of a turn each time until the valve cusps are destroyed.
5. If strong resistance is encountered, the valvulotome should be re-advanced, rotated a quarter of a turn and the pass attempt repeated.
6. The valves are deemed to be destroyed when pulsatile blood flow is observed distally. The distal anastomosis may then be performed.
7. Operative angiography is recommended for arterio-venous fistulas identification purposes and to confirm valve disruption and patency.

WARRANTY

Cardial warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Cardial's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL CARDIAL BE LIABLE TO YOU FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

References:

1. J.M. Chevalier et al. : Saphenous in situ bypass in critical limb ischaemia. Abstract from "Spinal cord simulation for peripheral vascular diseases. Advances and controversies". Libro del ano S.L: chapter 15, 1994.
2. J.M. Chevalier, B. Enon, J. Pillet : Destruction des valvules de la saphène interne dans la technique du pontage fémoro-jambier « in situ ». Nouv Press Méd 9 : 39-40, 1980.

TRADEMARKS

Bard is a registered trademark of C. R. Bard, Inc., or an affiliate.

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