**MollRing Cutter<sup>®</sup> Transection Device**

Instructions for Use - English

**MollRing Cutter<sup>®</sup> Transection Device**

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**MollRing Cutter Transection Device**

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**MollRing Cutter<sup>®</sup> Transection Device**

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**MollRing Cutter<sup>®</sup> Transection Device**

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**MollRing Cutter<sup>®</sup> Transection Device**

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**MollRing Cutter<sup>®</sup> Transection Device**

## MollRing Cutter® Transection Device

(Model Numbers 4200-40, 4200-41, 4200-42, 4200-43, 4200-44, 4200-45, e4200-40, e4200-41, e4200-42, e4200-43, e4200-44, and e4200-45)

Remote Endarterectomy Device

Instructions for Use - English

**Rx only** **STERILE EO** 

This device should be used by physicians trained in vascular surgery, specifically endarterectomy procedures.


CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

### Description

The LeMaitre® MollRing Cutter® Transection Device consists of two coaxial metal rings which are connected to a handle through a long rod and tube.

The device is designed to be introduced over the proximal end of the atheromatous core of a peripheral vessel and advanced to its distal end. The distal end of the atheromatous core is transected from the vessel wall by activating the cutting mechanism on the handle of the device. This selectively displaces the proximal and distal rings of the device so that they are not aligned with each other bringing the two cutting edges together transecting the core.

### Contents of Sterile Package

One MollRing Cutter Transection Device sterilized by ethylene oxide (EO). See package labeling for usable length and ring size.

### Specifications

See package labeling for product information.

Storage, shipping conditions: -10° to 50° C, up to 95% relative humidity.

Proper stock rotation should be employed.

### Indications For Use

This device is intended for remote intraluminal stripping of peripheral blood vessels during vascular reconstruction. General indications for peripheral vessel reconstruction include disabling claudication, ischemic pain at rest, or ischemic skin lesions or gangrene.

### Contraindications

Any evidence of aneurysmal change or severe calcification may contraindicate endarterectomy since continued aneurysmal degeneration of the adventitia may be expected.

In some cases, the vessel wall layers may be quite coherent and the plaque segment may be difficult to separate from the vessel wall. When the vessel is in this state of disease, it is often advisable to abandon endarterectomy and consider an alternative method of vascular reconstruction.

Use of this device is contraindicated in the coronary, internal carotid arteries, and other vessels of the central circulatory system.

### Precautions

The MollRing Cutter Transection Device should be advanced through the vessel and should cut the atheroma only under fluoroscopic guidance. Always verify the correct positioning of the rings.

Do not use excessive force with this device to strip and/or remove plaque segments from an occluded vessel.

This device is not designed to cut calcified material. Attempted transection of such material may damage the MollRing Cutter Transection Device.

If the rings, rods, or tubing of the MollRing Cutter Transection Device become bent, do not attempt to repair the device or use in the patient. Rather withdraw the device from the patient and dispose of the device.

As with any procedure, extreme care and discretion must be exercised at all times.

The physician must ensure that all excised atheromatous material is removed from the patient.

Complications associated with PTA and the use of endarterectomy devices include:

- Incisional/puncture site inflammation, hematoma or bleeding
- Wound infection (incisional)
- Vessel injury, dissection, perforation or rupture
- Vessel spasm
- Emboli, thrombosis or acute occlusion
- Fever, drug reactions or allergic reactions to contrast medium
- Conversion to bypass
- Limb loss due to distal embolization
- Renal failure
- Death

Devices are STERILE unless the package is opened, damaged or broken or the expiration date has passed.


### **Warnings**

This device is for single use only. Do not resterilize or reuse.

### **Instructions For Use**

The exact technique for endarterectomy is determined by the physician. One suggested technique is as follows:

1. Make an incision to expose the proximal aspect of the vessel to be treated.
2. Make an appropriate arteriotomy at the proximal end of the vessel to expose the atheroma.
3. Make a slanting, bevel-type incision of the proximal atheroma to start a cleavage plane between the atheroma and surrounding media.
4. Insert a conventional ring stripper into the vessel and advance it to a point distal to the occluded segment. A gentle rotating and thrusting motion is used to facilitate the development of the cleavage plane. Remove the ring stripper following completion of an adequate cleavage plane.

5. The point of endarterectomy termination will be determined from pre-operative arteriography. Fluoroscopy is used during the endarterectomy procedure. Advance the MollRing Cutter Transection Device to the desired end point.

*PRECAUTION: This device is not designed to cut calcified material. Attempted transection of such material may damage the MollRing Cutter Transection Device. If the rings, rods, or tubing of the MollRing Cutter Transection Device become bent, do not attempt to repair the device or use in the patient. Rather withdraw the device from the patient and dispose of the device.*

6. The atheroma is cut by rotating the lock lever counter-clockwise until it stops and then, while holding the blue knob stable, pulling back the white handle away from the cutting rings until it stops, bringing the cutting rings to their maximum separation.
7. Typically the atheroma core is grasped at the arteriotomy end using forceps. With the MollRing Cutter Transection Device knob slightly advanced, simultaneously pull the atheroma with the forceps and withdraw the MollRing Cutter Transection Device to remove the atheroma core.
8. Flush the vessel with saline to clear any loose, residual material. Angiography, angioscopy or intravascular ultrasound may be used to verify vessel patency following the procedure.
9. Secure any intimal flaps that may possibly dissect when flow is reestablished.
10. Close the arteriotomy.
11. Note: Physician experience and discretion will determine the appropriate drug regimen for each patient.

### **Resterilization/Re-use**

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

### **Safe Handling and Disposal**

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

### **Cleaning:**

1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

### **Packaging:**

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.

4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803

**Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

## MollRing Cutter® Transection Device

(Modellnummern 4200-40, 4200-41, 4200-42, 4200-43, 4200-44, 4200-45, e4200-40, e4200-41, e4200-42, e4200-43, e4200-44, und e4200-45)

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LeMaitre and MollRing Cutter are registered trademarks
©2019 LeMaitre Vascular, Inc. All rights reserved.
LeMaitre Vascular GmbH
of LeMaitre Vascular, Inc.