Syntel® Silicone Embolectomy Catheter - Spring Tip

Instructions for Use - English

Syntel® Embolektomiekatheter aus Silikon - Federspitze

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Product Description

Syntel Embolectomy Catheters consist of a catheter body with a silicone balloon on the distal end and an inflation hub on the proximal end.

Indications for Use

Syntel Embolectomy Catheters are indicated for the removal of emboli and thrombi from the peripheral arterial system.

Contraindications

Syntel Embolectomy Catheters are contraindicated for endarterectomy procedures or for use in the venous system, grafts, shunts, or as a vessel dilator.

Instructions for Use

1. The catheter balloon should be inflated with a sterile, nonparticulate, blood compatible liquid or CO₂.
2. If using CO₂, corrections to inflation volume may be necessary due to CO₂ diffusion through the balloon wall.
3. Purge the catheter of all air with a sterile medium.
NOTE: Inspection for leaks should be performed during this step. Do not use if leakage is observed.
4. For optimal control, use the recommended syringe capable of holding inflation medium. See Specification Table.
5. Insert the catheter through the arteriotomy and position the balloon just past the occlusion.
6. Inflate the balloon with a sterile medium. Inflation should continue until the balloon can be felt to have engaged the vessel wall.
7. Withdraw the catheter and thrombus through the arteriotomy.

Warnings

- Air should not be used to inflate the balloon.
- Do not exceed maximum pull force specifications as balloon rupture and catheter separation may result (See Specification Table).
- Do not use if product or sterile packaging is damaged.

Potential Complications

Complications associated with the use of embolectomy catheters include, but are not limited to: systemic infection, local hematomas, intimal disruptions, arterial dissection, perforation and rupture, hemorrhage, arterial thrombosis, distal emboli of blood clots or arteriosclerotic plaque, air embolus, aneurysms, arterial spasms, arteriovenous fistula formations, and balloon rupture or tip separation with fragmentation and distal embolization.

Precautions

- Balloon degradation and rupture may result from exposure to adverse environmental conditions, excessive handling, or deposits within the vessel.
- Maximum recommended inflation volumes should not be exceeded (See Specification Table).
- Balloon rupture is sensed by a decrease in resistance of the syringe plunger during inflation. If a balloon ruptures, discontinue inflation and remove the catheter at once.
- To minimize inner wall pressure and shear forces to the inner surface of the artery, use of the smallest inflated balloon diameter to remove the obstructing material is recommended.

Storage

Handle with care. Product should be stored in a clean, cool, dry area away from chemical fumes.

How Supplied

This device is supplied sterile and non-pyrogenic.

Restoration/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only.

Safe Handling and Disposal

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
 - a) Sodium hypochlorite solution (500-600 mg/l), or
 - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
 - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
 - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.

5. Shipping paper and content marking for the outer shipping container are not required.

6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803

Warranty Coverage

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Other than as explicitly provided herein, LEMAITRE VASCULAR MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR (AS USED IN THIS PARAGRAPH, INCLUDING ITS EMPLOYEES, OFFICERS, DIRECTORS, OR AFFILIATES) BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and Syntel are registered trademarks of LeMaitre Vascular, Inc.
© 2021 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH