**LeMills Valvulotome**

Instructions for Use - English

**LeMills Valvulotom**

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**Valvulotome LeMills**

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**Valvulotomo LeMills**

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**Valvulotomo LeMills**

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**Valvulótomo LeMills**

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## STERILE EO Rx only

### Description

The LeMills family of Valvulotomes are devices that cut venous valves during vascular procedures. The devices consist of stainless steel blades and plastic handles.

### Indication for Use

The LeMills Valvulotomes are intended to cut venous valves.

### Contraindications

1. endarterectomy procedures,
2. thrombolysis procedures,
3. vein stripping procedures, or
4. embolectomy procedures.

### Warnings

1. Do not use if the package or device is damaged.
2. Do not use if inner packaging is opened outside a sterile environment.
3. The LeMills Valvulotomes are to be used by a qualified physician.
4. Do not pass the LeMills Valvulotome through a vessel that has undergone synthetic grafting or contains implants.
5. Do not use device for valvulotomy unless the target vein is fully distended by arterial blood flow (due to anastomosis of artery to vein) or by saline injection.

### Precautions


1. United States Federal and other law restricts this device to sale on or by the order of a physician.
2. Do not use devices past the expiration date printed on the labeling.
3. These devices are single use only. Do not reuse, reprocess, or re-sterilize (reference "Re-sterilization/Repackaging" section of this document for further details). The device is NOT autoclavable.
4. Keep any dissected portion of the vein straight and avoid twisting the vein during valvulotomy.
5. Use caution when tracking the LeMills Valvulotomes in thrombophlebitic veins.

### Potential Complications


- Vessel wall perforation
- Intimal disruption
- Entrapment of cutting blades in branches or anastomoses
- Inadequate valve disruption
- Vein Stricture
- Thrombus formation
- Post-procedure wound necrosis

### Specifications

Antegrade LeMills Valvulotome:

Retrograde LeMills Valvulotome:

Catalog Number	Description	Useable Length (cm)
1050-00, e1050-00	Retrograde LeMills Valvulotome	12.5
1050-01, e1050-01	Antegrade LeMills Valvulotome	12.5
1050-02, e1050-02	Antegrade & Retrograde LeMills Valvulotomes	12.5; 12.5

### To Open Package

1. Open box and remove sealed sterile unit.
2. Open sealed sterile unit and present contents to personnel in sterile environment.

### Pre-Use Check

Inspect blades for damage.

### Instructions for Use

1. Insert the device into the vein to the treatment area.
2. If using Retrograde LeMills, slowly retract the device to disrupt the valve sets. If using Antegrade LeMills, slowly advance the device to disrupt the valve sets.
3. Steps 1 and 2 may be repeated as necessary.
4. Confirm the free flow of blood through the vein.

### Storage

Store in a cool, dry place.

### Reste r ilization/Re-use

These devices are single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

### Safe Handling and Disposal

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed accorded to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no heavy metals or radioisotopes, and is not infectious or pathogenic. Please consult local regulations to verify proper disposal of sharps.

### Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
  - (1) Sodium hypochlorite solution (500-600 mg/l), or
  - (2) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - (1) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - (2) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

### Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be

- labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
  4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
  5. Shipping paper and content marking for the outer shipping container are not required.
  6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803

### **Limited Product Warranty; Limitation of Remedies**

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre is a registered trademark of LeMaitre Vascular, Inc.
©2020 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH