**Pruitt<sup>®</sup> Aortic Occlusion Catheter**

English — Instructions for Use

**Pruitt<sup>®</sup> Aortenokklusionskatheter**

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**Pruitt® Aortic Occlusion Catheter (2100-12)**  
English – Instructions for Use

**STERILE** |  **Rx only**

**Introduction**

Gaining control of bleeding of the proximal abdominal aorta with an aortic aneurysm is a difficult procedure requiring rapid occlusion of the vessel. LeMaitre Vascular's Aortic Occlusion Catheter, designed and tested by surgeons, offers the option to accomplish this via a transluminal approach when there is an acute aortic emergency, or when dissection and cross-clamping are difficult and time consuming (i.e. when the aneurysm is excessively large, dense, adherent to neighboring structures, such as the intestines, or especially if it is ruptured.) Intraluminal balloon occlusion may be accomplished by direct insertion through the wall of the aneurysm.

**Product Description**

The Aortic Occlusion Catheters are 12 French (4.0 mm), dual lumen catheters with a large, latex balloon (maximum liquid inflation capacity 50 mL) specifically designed and sized for use in the outlined general procedures. The first lumen (inflation lumen indicated by the white stopcock) is used for balloon inflation, while the second lumen (irrigation lumen indicated by the blue stopcock) allows access to the vessel distal to the occlusion. Other features include 2 stopcocks with a luer-lock fitting at the proximal end of the irrigation lumen to facilitate control of such procedures, a balloon wall thickness designed to reduce the possibility of puncture by calcium deposits, and a stopcock to maintain balloon inflation level throughout the procedure.

A stainless steel stylet is inserted in the irrigation lumen of the catheter and serves as a stiffening medium to aid the physician during introduction of the catheter into the patient's aorta.

**How Supplied**

The Aortic Occlusion Catheter with stylet, Model 2100-12 is sterile and nonpyrogenic in sealed peel-open packages. Sterility is assured as long as the packages are unopened and undamaged.

**Indication**

The Pruitt Aortic Occlusion Catheter and kit are designed for the purpose of obtaining rapid control of in-flow blood in the abdominal aorta in cases of ruptured aortic aneurysm or in other conditions when dissection of the neck of the aneurysm for different reasons may be especially difficult. This application of direct internal balloon occlusion is in lieu of external, proximal cross-clamping of the abdominal aorta and may prevent technical complications from special physiological situations. Intraluminal balloon occlusion may be accomplished by direct insertion through the wall of the aneurysm.

**Contraindications**

1. The catheter is not to be used as a dilation catheter.
2. The catheter is not to be used for the introduction of drugs other than saline, heparin, and contrast media.
3. The catheter is a temporary device and cannot be implanted.

**Warnings**

1. Do not reuse. The catheter is for single use only.
2. Air or gas should not be used to inflate the balloon during patient use.
3. Do not inflate the balloon to any greater volume than is necessary to obstruct the blood flow. DO NOT EXCEED the recommended maximum balloon inflation capacity (maximum liquid inflation capacity 50 mL).
4. Exercise caution when encountering extremely diseased vessels. Arterial rupture or balloon failure due to sharp calcified plaque, may occur.
5. Deflate the balloon prior to inserting or withdrawing the catheter. Avoid using excessive force to push or pull catheter against resistance.
6. The possibility of balloon rupture or failure must be taken into account when considering the risk involved in a balloon catheterization procedure.
7. All agents to be infused should be used according to the manufacturer's Instructions for Use.

**Precautions**

1. Inspect the product and package prior to use and do not use the catheter if there is any evidence that the package or the catheter has been damaged.

2. Avoid extended or excessive exposure to fluorescent light, heat, sunlight, or chemical fumes to reduce balloon degradation. Excessive handling during insertion, or plaque and other deposits within the blood vessel may damage the balloon and can increase the possibility of balloon rupture.
3. Ensure proper connections between all syringes and hubs to avoid the introduction of air.
4. Do not grasp the balloon with instruments at any time to avoid damage to the latex.
5. Aspirate the irrigation lumen of the catheter during insertion until there is free back flow of blood from the catheter to reduce the chance of air embolism.

**Adverse Events**

As with all catheterization and surgical procedures, complications may occur. These may include, but are not limited to:

- infection
- local hematomas
- intimal disruption
- arterial dissection
- vessel perforation and rupture
- hemorrhage
- arterial thrombosis
- distal emboli of blood clots or arteriosclerotic plaque
- air embolus
- aneurysms
- arterial spasm
- arteriovenous fistula formation
- balloon rupture
- tip separation with fragmentation and distal embolization
- renal insufficiency
- paraplegia

**Procedure**

Pretest (perform before patient use)


1. Aspirate the balloon completely prior to inflation of the balloon with liquid.
2. Inflate the balloon with sterile saline and inspect for leaks. If there is any evidence of leaks around the balloon or if the balloon will not remain inflated, do not use the product.
3. Check the balloon by inflating and deflating with sterile saline for injection before use. If the balloon does not appear to function normally, do not use the product.

**General Use**

1. Approach the abdominal aorta through a mid-line laparotomy or through any other approach of the surgeon's choice. Depending on the urgency dictated by the situation, dissection of the aortic bifurcation or distal iliac arteries may be carried out.
2. Insert an 18-gauge needle into the center of the aneurysm (Figure A). Blood should pulse from the needle hub indicating that the needle is indeed located in the bloodstream. A 5 mL syringe may be used to aspirate blood flow or clear any thrombus that may occlude the needle during its placement into the true lumen of the aorta.
3. Pass a "J" guidewire (0.9mm x 55 cm) through the needle and maneuver it through the neck of the aneurysm into the suprarenal aorta (Figure B).
4. Remove the needle and insert the introducer sheath (7.3 mm) over the guidewire into the aortic lumen (Figure C).

NOTE: The balloon catheter may be left in loco during graft placement and the proximal anastomosis carried out around the balloon catheter. If such a choice is made, thread the catheter through the graft BEFORE insertion (Figure E and F).

5. Remove the dilator from the introducer sheath and thread the balloon catheter over the "J" guidewire up through the outer sheath of the introducer (Figure D).
6. Flush irrigation lumen (indicated by the blue stop-cock) with saline for irrigation to remove air.
7. Insert the balloon catheter into the suprarenal aorta, remove the guidewire and inflate the balloon. To temporarily occlude the aorta during surgery, position the catheter balloon (white stopcock indicates the balloon inflation lumen) with the aortic lumen to the point requiring the occlusion. Inflate the balloon with sterile saline for injection (or other sterile, blood compatible fluid that may include highly dilute and nonparticulate radiopaque solutions) using a 30 mL syringe to occlude the vessel, taking care not to over inflate the balloon. Do not inflate the balloon to any greater volume than necessary to obstruct the blood lumen. Close the inflation stopcock to maintain balloon inflation. During positioning, aspirate the irrigation lumen until there is free back flow of blood from the catheter, to reduce the chance of air embolism (1).
8. Once the catheter is suitably positioned, introduction or withdrawal of fluids to areas distal to the point of occlusion may be accomplished through the irrigation-lumen. This is facilitated through the luer-lock fitting at the base of the irrigation stopcock.
9. After the proximal anastomosis is completed, withdraw the catheter, clamp the proximal end of the graft and complete the distal graft anastomoses (Figure G and H).  
or  
Perform the dissection of the neck of the aneurysm under the protection of the balloon catheter, remove the balloon catheter and clamp the neck itself (Figure I and J). The catheter is occluding blood flow to the kidneys; therefore, it should not be left in longer than 30-45 minutes.

#### Storage Shelf Life

The shelf life is indicated by the USE BY date on the package label. The Use By date printed on each label is NOT a sterility date. The Use By date is based on the normal life expectancy of the natural latex balloon when properly stored. The use of the catheter beyond the expiration date is not recommended because of potential balloon deterioration. LeMaitre Vascular, Inc. does not make provisions for replacing or reprocessing expired product.

Since natural rubber latex is affected by environmental conditions, proper storage procedures must be practiced to achieve optimum shelf life. The product should be stored in a cool dark area away from fluorescent lights, sunlight, and chemical fumes to prevent pre-mature deterioration of the rubber balloon. Proper stock rotation should be practiced.

#### Re-sterilization/Repackaging

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

#### Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME.

LeMaitre Vascular makes no representation regarding the suitability for any particular treatment in which this device is used, which determination is the sole responsibility of the purchaser. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

#### References

1. Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
2. Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978: 83:549.
3. Dujovny, M, Laha RK, Barrionuevo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7:39.
4. McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
5. Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Jour-nal of Cardiovascular Sur-gery* 1987: 28:549-551.
6. Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
7. Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Literaturhinweise

1. Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. Surgery. 1980;634-638.
2. Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter, Surgery. 1978: 83:549.
3. Dujovny, M, Laha RK, Barrioneuo, P. Endothelial changes secondary to use of the Fogarty Catheter. Surg, Neuro, 1977: 7:39.
4. McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
5. Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. Jour-nal of Cardiovascular Sur-gery 1987: 28:549-551.
6. Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. Surgery. 1970: 68(5):774-777
7. Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. Surgical Rounds. 1988:17-27.

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#### Références

1. Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
2. Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978: 83:549.
3. Dujovny, M, Laha RK, Barrionuevo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7:39.
4. McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
5. Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Journal of Cardiovascular Sur-gery* 1987: 28:549-551.
6. Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
7. Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Bibliografia

1. Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
2. Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978; 83:549.
3. Dujovny, M, Laha RK, Barrionero, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977; 7:39.
4. McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
5. Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Journal of Cardiovascular Sur-gery* 1987; 28:549-551.
6. Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
7. Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Referencias

- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978; 83:549.
- Dujovny, M, Laha RK, Barrionuevo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7:39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Journal of Cardiovascular Sur-gery* 1987; 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
- Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Bibliografia

1. Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;63:4-638.
2. Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978; 83:549.
3. Dujovny, M, Laha RK, Barrioneu, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7-39.
4. McCughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
5. Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Jour-nal of Cardiovascular Sur-gery* 1987: 28:549-551.
6. Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
7. Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Referanser

1. Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. Surgery. 1980;634-638.
2. Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter, Surgery. 1978: 83:549.
3. Dujovny, M, Laha RK, Barrionuo, P. Endothelial changes secondary to use of the Fogarty Catheter. Surg, Neuro, 1977: 7:39.
4. McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
5. Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. Jour-nal of Cardiovascular Sur-gery 1987: 28:549-551.
6. Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. Surgery. 1970: 68(5):774-777
7. Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. Surgical Rounds. 1988:17-27.

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#### Referencer

- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978; 83:549.
- Dujovny, M, Laha RK, Barrionuevo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7:39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Jour-nal of Cardiovascular Sur-gery* 1987: 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
- Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Referenser

- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978: 83:549.
- Dujovny, M, Laha RK, Barrioneuo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7:39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Jour-nal of Cardiovascular Sur-gery* 1987: 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
- Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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#### Viitteet

- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;634-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. *Surgery*. 1978: 83:549.
- Dujovny, M, Laha RK. Barrioneuo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977: 7:39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Jour-nal of Cardiovascular Sur-gery* 1987: 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
- Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. Surgery. 1980;634-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter. Surgery. 1978; 83:549.
- Dujovny, M, Laha RK. Barrionero, P. Endothelial changes secondary to use of the Fogarty Catheter. Surg, Neuro, 1977; 7:39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. Jour-nal of Cardiovascular Sur-gery 1987; 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. Surgery. 1970: 68(5):774-777
- Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. Surgical Rounds. 1988:17-27.

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- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. Surgery. 1980;634-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter, Surgery. 1978: 83-549.
- Dujovny, M, Laha RK. Barriouneo, P. Endothelial changes secondary to use of the Fogarty Catheter. Surg, Neuro, 1977: 7-39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. Jour-nal of Cardiovascular Sur-gery 1987: 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. Surgery. 1970: 68(5):774-777
- KunkeI J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. Surgical Rounds. 1988:17-27.

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- Burdick JF, Williams GM. A study of the lateral wall pressure exerted by balloon-tipped catheter. *Surgery*. 1980;63:4-638.
- Chidi D, DePalma, RG. Atherogenic potential of the embolectomy catheter, *Surgery*. 1978; 83:549.
- Dujovny, M, Laha RK. Barrioneuo, P. Endothelial changes secondary to use of the Fogarty Catheter. *Surg, Neuro*, 1977; 7:39.
- McCaughan, JJ, Young, JM. IntraArterial Occlusion in Vascular Surgery. 1970: 695-703. Rev. A
- Robicsek F, Pruitt JC. Transaneurysmal aortic balloon tamponade. *Jour-nal of Cardiovascular Sur-gery* 1987: 28:549-551.
- Robicsek F, Daugherty HK, Mullen DC. The elective use of balloon obstruction in aortic surgery. *Surgery*. 1970: 68(5):774-777
- Kunkel J, Moore WS. Man-agement Errors in the Treatment of Ruptured Aortic Aneurysm. *Surgical Rounds*. 1988:17-27.

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LeMaitre Vascular, Inc.  
63 Second Avenue  
Burlington, MA 01803

LeMaitre and Pruitt are registered trademarks of LeMaitre Vascular, Inc.

© 2010 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.

R1577-00 Rev. E 01/10

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LeMaitre Vascular GmbH