Omniflow® II Vascular Prosthesis

Instructions for Use - English

Omniflow® II Gefäßprothese

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Omniflow® II Protesi vascolare

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Omniflow® II Prótesis vascular

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Omniflow® II Prótese vascular

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Description

Omniflow® II Vascular Prosthesis is constructed of polyester mesh within an ovine fibrocollagenous tissue matrix. The prosthesis is sterilised using aseptic processing techniques. The prosthesis is also treated with a glutaraldehyde process which crosslinks the collagen fibers and minimizes antigenicity.

The prosthesis is supplied sterile and nonpyrogenic in a solution of 50% ethanol. The prosthesis remains sterile unless the primary package is opened or damaged. -

The Omniflow II Straight Vascular Prosthesis is mounted on a plastic mandrel contained in a plastic tube. The mandrel design prevents the prosthesis from slipping off the mandrel when it is removed from the tube. The diameter and minimum length of the prosthesis is specified on the label applied to the tube. The plastic tube is sealed in a non-sterile outer plastic bag.

Note: The graduations on the Omniflow tube are for reference only. The tube is not to be used as a measuring device.

The Omniflow II Curved Vascular Prosthesis is contained in a sterile flexible inner pouch within a peelable foil outer pouch. The diameter and minimum length of the prosthesis is specified on the label applied to the the outer surface of the outer pouch.

The Omniflow II Vascular Prosthesis is considered MR safe.

Indications for Use

The Omniflow II vascular prosthesis is indicated for the replacement, reconstruction, bypassing or patching of diseased vessels in patients suffering occlusive or aneurismal diseases, in trauma patients requiring vascular replacement or for patients requiring vascular access such as for haemodialysis.

The Omniflow II Straight Vascular Prosthesis is intended:

1. To bypass, replace or reconstruct diseased or injured blood vessels
2. To patch and repair peripheral vessels
3. For arteriovenous access when a straight configuration is required.

The Omniflow II Curved Vascular Prosthesis is intended for arteriovenous access when a looped configuration is required.

Warnings

1. Do NOT re-sterilise the Omniflow II prosthesis. Use the prosthesis immediately after opening the package and discard any unused portions.
2. Do NOT use the prosthesis if the primary package is damaged as sterility may be compromised.
3. Do NOT use the prosthesis if it is not completely covered by the storage solution.
4. Do NOT use the prosthesis if the solution level in the vertical position is below the anti-roll nubs on the tube.
5. Do NOT attempt to reposition the prosthesis after removal of the tunnelling instrument.
6. Do NOT straighten the curved prosthesis during preparation or implantation, as this will cause disruption of the mesh tissue interface.
7. Do NOT use the straight prosthesis to fashion a looped arteriovenous access as this may cause kinking.
8. Do NOT pull, stretch, twist, squeeze or pinch the body of the prosthesis.
9. Do NOT use ablation techniques such as cutting balloons, laser, or radio frequency ablation with the Omniflow II prosthesis.
10. Do NOT attempt to dilate the prosthesis with balloon angioplasty or stenting procedures.
11. The Omniflow II prosthesis should only be implanted by trained surgeons.
12. The use of the Omniflow II prosthesis in the coronary artery has not been evaluated.
13. Ethanol is a highly flammable liquid and vapor. Keep away from heat, sparks, and open flames.

Technical Information/Precautions

1. Ensure the rinsing procedure has been performed to remove the storage solution prior to implanting the prosthesis. Failure to do so may cause occlusion. Keep the prosthesis moist with sterile physiological saline during the procedure.
2. The use of a hollow tunnelling instrument for the passage of the prosthesis is essential. Failure to do so may cause disruption to the biosynthetic material and lead to occlusion, dilatation or aneurysm formation. The inner diameter of the tunneler should be at least 3mm larger than the indicated inner diameter of the prosthesis.
3. Ensure that the prosthesis does not become twisted when passing through the tunnelling instrument as this may lead to occlusion.
4. Avoid cross clamping with metal instruments as this may damage the prosthesis and cause occlusion, dilatation or aneurysm formation. If clamping is necessary use only atraumatic clamps and avoid repeated or excessive clamping in the same position on the prosthesis.
5. The prosthesis has minimal longitudinal elasticity. Ensure the prosthesis is cut to the correct length. If it is too short it may cause suture pullout with a risk of anastomotic aneurysm. If it is too long it may kink and cause occlusion.

6. Cut off the sections of the prosthesis which were clamped during rinsing. Ensure that the full wall thickness and a mesh eyelet are incorporated with each stitch when performing the anastomosis. Failure to do so may result in stitch pullout and anastomotic aneurysm formation.
7. Do not implant Omniflow II into a site with an active infection unless the surgeon determines there is not a more suitable alternative for preventing amputation or death.
8. When the prosthesis is used for arteriovenous access some redness and swelling may be present over the implant area for a few days following implant.
9. Insufficient data is available on which to base any conclusions regarding the use of the Omniflow II vascular prosthesis for aortocoronary bypass procedures.
10. Failure to heparinize the prosthesis (i.e. in the case of patients who cannot tolerate heparin) may result in a higher likelihood of thrombosis or occlusion post-implantation, the extent of which has not been established.
11. Omniflow II cannot grow in diameter or length and thus should not be implanted into infants or children unless a plan for its replacement is established and that no other suitable alternative treatment option exists.

Thrombectomy

Prompt intervention when occlusion occurs has been found to restore prosthetic function. Therefore the patient should be advised to report to the clinician if symptoms return. When performing a thrombectomy choose an embolectomy catheter of the appropriate size and follow the catheter manufacturer's instructions for use.

To avoid damage to the wall and flow surface of the prosthesis:

1. Do NOT exert undue pressure when inflating the balloon
2. Do NOT exert undue force when withdrawing the catheter from the prosthesis.

Contraindications

The prosthesis should not be used in patients with a known hypersensitivity to ovine material or glutaraldehyde.

Potential Complications with the use of Vascular Prostheses

Complications may occur with the use or in conjunction with any vascular prosthesis and include but are not limited to: infection, thrombosis/ occlusion, dilatation, leakage, aneurysm formation, pseudoaneurysm formation, suture pullout and adverse tissue responses.

Late aneurysm formation in Omniflow II has been reported (more than 4 years after implantation). Ongoing patient monitoring is recommended.

Wall integrity of the prosthesis may be adversely affected by collagenase-producing microorganisms.

Pre-implant preparation

Prepare a sterile basin with adequate sterile physiological saline to cover the prosthesis, a sterile 20 ml syringe, sterile heparin 5,000 I.U. per mL.


Removing the prostheses from the container

STRAIGHT VASCULAR PROSTHESIS

1. Remove the bag and tamper seal from the tube. The outer bag is not a sterile barrier and the exterior of the tube, cap and stopper are non-sterile.
2. Check the storage solution level of the tube in the vertical position. The solution level must be above the top of the graft and above the bottom of the anti-roll nubs in order to ensure proper prosthesis hydration.
3. Remove the cap and discard. Remove the stopper with forceps and discard. The forceps must also be discarded as they are now considered non-sterile. Fig. 1.
4. The mandrel will float above the opening of the tube. Aseptically remove the prosthesis from the tube by grasping the end of the mandrel and lifting it gently out of the tube. Fig. 2. Do not allow the prosthesis to touch the outer top edge of the tube. The "T" on the bottom of the mandrel will prevent the graft from slipping off the mandrel.
5. Grasp the mandrel by the T, turn it upside down and allow the prosthesis to slide into the basin of sterile saline. If the prosthesis does not slide off easily, gently push it down the mandrel. Do NOT pull the prosthesis. Discard the mandrel and tube. Discard the ethanol storage solution according to local guidelines.

CURVED VASCULAR PROSTHESIS

1. Aseptically remove the sterile inner pouch containing the prosthesis from the outer foil pouch. Fig. 3. Only the inner pouch is sterile.
2. Cut the corner of the inner pouch with sterile scissors and allow the storage solution to drain into a container. Fig. 4. Discard the container of ethanol storage solution according to local guidelines.
3. Remove the prosthesis from the inner pouch and place it into the basin of sterile physiological saline.

Rinsing the prosthesis


STRAIGHT VASCULAR PROSTHESIS

1. Use a sterile syringe to flush the lumen of the prosthesis with 20 ml of sterile physiological saline. Fig. 5. Repeat at least twice.
2. Clamp one end of the prosthesis. Insert sterile heparin at a concentration of 5,000 I.U. per mL into the lumen of the prosthesis so that the entire inside surface of the graft is coated with the solution (50,000 I.U. in 10 mL is sufficient for any length prostheses).
3. Apply a clamp to the other end of the prosthesis.
4. Place the prosthesis back into the basin completely covered by the sterile saline until it is time for implantation. This should be a minimum of 10 minutes.

When it is time to implant the prosthesis, remove the clamps from the prosthesis and allow the heparin and saline to drain out.

CURVED VASCULAR PROSTHESIS

The curved prosthesis is rinsed in the same manner as described above for the straight prosthesis. The curved configuration is maintained throughout the procedure. Fig. 6. Do NOT straighten the prosthesis.

Implantation

Prepare the incisions in the normal fashion.

STRAIGHT VASCULAR PROSTHESIS

Form a tunnel for the prosthesis using a metal or plastic hollow tunnelling instrument and pass the prosthesis through the instrument. Fig. 7.

CURVED VASCULAR PROSTHESIS

Marking the proposed graft layout on the skin with a sterile surgical marker may facilitate final positioning of the prosthesis. Maintain the curve of the graft by inserting 2 hollow tunnelling instruments and passing both straight sections of the graft through simultaneously from the distal to the proximal incision. Fig. 8. Ensure that the graft does not kink at the curve during placement.

Carefully remove the tunnelling instrument when the graft has been correctly placed.

Anastomoses

1. When fashioning the ends for the anastomoses it is essential to cut off those portions of the prosthesis that were clamped, as clamping can damage the tissue matrix. No special cutting technique is required given the biosynthetic material of the graft..
2. Monofilament polypropylene is the preferred suture material. During suturing evert the edges of the prosthesis and ensure the full wall thickness and a mesh eyelet are taken up with each stitch. Fig. 9. It is advised to only use non-cutting needles when suturing.
3. Prior to completion of the second anastomosis flush the prosthesis with blood to ensure complete removal of the heparin.
4. An intra-operative angiogram may be performed at the time of the procedure to document function. Ensure that the artery rather than the prosthesis is used for the injection.

At the completion of the procedure the patient record labels supplied with the prosthesis should be affixed in the patient and hospital records.


Fig. 9

Arteriovenous Access

1. It is preferable to allow the prosthesis to heal for a minimum of two weeks prior to commencement of puncturing.
2. Do NOT puncture the arteriovenous access repeatedly in the same site as this may lead to disruption of the wall of the prosthesis, haematoma formation or pseudoaneurysm formation. Rotation of the puncture sites along the length of the prosthesis is necessary.
3. Do NOT puncture the curve of the prosthesis.
4. Do NOT puncture within 3 cm of the anastomotic sites.
5. Strict adherence to aseptic technique is essential to minimise the risk of infection.
6. Moderate compression of the puncture site will achieve haemostasis.

Vascular Patch

Cut the tubular prosthesis longitudinally to fashion into a vascular patch. Ensure that the blood flow surface is not damaged. When suturing the patch into place, ensure that the full wall thickness and the mesh are incorporated into each stitch.

Patient Information

The biosynthetic nature of Omniflow II does not change advice given to patients regarding post-operative care of the prosthesis as compared to other typical vascular prostheses. As such, good practice should include the following:

For bypass patients, avoid any compression of the prosthesis either by overly tight wrapping or clothing or laying on the prosthesis.

Patients should report any abnormal signs immediately to the treating surgeon to reduce the risk of complications.

An Omniflow Patient Implant Card is provided with each prosthesis. Please complete the front side of the card and provide it to the patient.

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1. Check the prosthesis every day. Feel for the pulse (thrill) in the prosthesis.
2. Inspect daily for signs of infection e.g. swelling, redness and pus.
3. Inspect daily for bruising or pseudoaneurysm formation.
4. Keep the skin over the prosthesis clean to help prevent infection.
5. Do not wear tight clothing or jewelry or drape or carry heavy items over the prosthesis as these may restrict blood flow.
6. Do not sleep on the arm with the prosthesis as this may occlude the blood flow.
7. Do not allow blood pressure testing or use of intravenous drips or the taking of blood from the arm with the prosthesis.
8. Patients should report any abnormal signs immediately to the dialysis unit to reduce the risk of complications.

Storage

The prosthesis must be kept in its original packaging until used. It must be kept in a dust free and dry environment at room temperature.

Safe Handling and Disposal

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located. This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal. Dispose storage solution according to local and federal regulations. Solution should not be disposed using septic systems. If there are no disposal restrictions, solution may be diluted and disposed in a sanitary sewer system. For more information see:

<https://www.osha.gov/Publications/glutaraldehyde.pdf>

Packaging and shipping of explanted Omniflow:

Return of the shipment to LeMaitre Vascular depends on 3 crucial questions:

1. Is the explant obtained from a patient with a known or presumed pathogenic condition at the time of the explant?
2. Is the explant obtained from a patient that has a known treatment history that involves therapeutic radionuclides within the last 6 months?
3. Has the clinician obtained consent from the patient for the sample to be returned to the manufacturer for research purposes?

In the event that the response to question 1 or 2 is affirmative, LeMaitre Vascular does not provide adequate guidance for shipment. THESE EXPLANTS SHOULD NOT BE RETURNED TO LEMAITRE VASCULAR UNDER ANY CIRCUMSTANCES. In these cases, the explant should be disposed according to local regulations.

For explants that do not have pathogenic or radiologic hazards, please use the following:

Pre-explantation:

1. If possible, perform a CT or Ultrasound scan of the device to document patency.
2. LeMaitre Vascular can accept clinical information that is patient-anonymized. LeMaitre Vascular requests information including:
 - a) The original diagnosis which resulted in use of the implant.
 - b) The patient's medical history relevant to the implant, including the hospital or clinic at which the device was implanted.
 - c) The patient's implant experience prior to implant removal.
 - d) The hospital or clinic at which the explantation was done, and date of retrieval.

Explantation:

1. Explanted Omniflow prostheses should be directly transferred to a sealable container filled with a solution of alkaline buffered 2% glutaraldehyde or 4% formaldehyde prior to shipment.
2. Cleaning of explanted prostheses should be minimal if necessary. Proteolytic digestion should not be used in any circumstances.
3. Omniflow explants should not be decontaminated under any circumstances. DO NOT autoclave the sample or use ethylene oxide gas to decontaminate.

Packaging:

1. Explants should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. Material that is absorbent and is cushioning should be selected for isolating the sealable container inside the secondary packaging. Primary and secondary packaging must then be packaged inside an outer package.
2. Explants in sealed primary containers should be labeled with an ISO 7000-0659 Biohazard symbol. The same symbol should be attached to the secondary packaging and to the outer packaging. Outer packaging should also be labelled with Name, Address and Telephone Number of Sender, and the statement, "Upon discovery of damage or leakage, the package should be isolated and sender notified".
3. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803, USA

Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. The unused section of the Omniflow II prosthesis or the explanted prosthesis should be considered as biohazard and be disposed per hospital policy or local regulations.

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and Omniflow are registered trademarks of LeMaitre Vascular, Inc.
©2022 LeMaitre Vascular, Inc. All rights reserved.
LeMaitre Vascular GmbH
LeMaitre Vascular, Ltd.