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Syntel™ silikoninen Over-the-Wire -embolektomiakatetri

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Syntel™ Silicone Over-the-Wire Embolectomy Catheter

Syntel™ Silicone Over-the-Wire Embolectomy Catheter

(Model Numbers A4E01, A4E02, A4E03, A4E04, A4E05, A4E06, A4E08, A4E09)

Instructions For Use – English

STERILE EO  **Rx only**

Product Description

The Syntel™ Silicone Over-the-Wire Embolectomy catheter consists of a double-lumen catheter body with a silicone balloon on the distal end and two standard luer connectors on the proximal end. One lumen is used for balloon inflation and the other is used for infusion of fluids, or insertion over a guidewire.

The Syntel™ Silicone Over-the-Wire Embolectomy catheter is marked every 10cm so that intraluminal length of the inserted portion of the catheter can be identified. Balloon inflation capacity is printed on the proximal end of the catheter. A syringe is included for balloon inflation. Marker bands are located at the proximal end and at the distal end of the balloon to aid in identifying balloon placement under fluoroscopy.

Indication for Use

The Syntel™ Silicone Over-the-Wire Embolectomy catheters are indicated for removal of thromboemboli from the peripheral arterial system, for occlusion of the vessel, and for infusion of fluids into a vessel.

Contraindications

The Syntel™ Silicone Over-the-Wire Embolectomy catheters are contraindicated for endarterectomy procedures or for use in the venous system.

Warnings

1. Do not use if the product or sterile package is damaged.
2. Do not reuse. The catheter is for single use only.
3. Do not use air to inflate the balloon during patient use.
4. Do not inflate the balloon to any greater volume than is necessary to obstruct the blood flow. DO NOT EXCEED the recommended maximum balloon inflation capacity (see specifications).
5. Exercise caution when encountering extremely diseased vessels. Arterial rupture or balloon failure due to sharp calcified plaque may occur. The possibility of balloon rupture must be taken into account when considering the risks involved in the catheterization procedure.
6. Deflate the balloon prior to catheter insertion. Avoid using excessive force to push or pull catheter against resistance.
7. Use all infusion agents according to their manufacturer's Instructions for Use.

Precautions

1. Maximum recommended inflation volumes should not be exceeded (See Specification Table).
2. Do not grasp the balloon with instruments at any time to avoid damage to the balloon.

Procedure

Pretest (Perform before patient use)

1. Inspect each package and do not use if there is any evidence that the package has been opened or damaged.
2. Aspirate the balloon completely prior to inflation.
3. Inflate the balloon with sterile saline and inspect for leaks. If there is any evidence of leaks around the balloon, if the balloon will not remain inflated or if the balloon does not appear to function normally, do not use the product.

CAUTION: Do not exceed the maximum recommended volume, as over inflation increases the possibility of balloon rupture. Refer to specifications for maximum inflation capacities.

General Use

1. The catheter balloon should be inflated with a sterile, nonparticulate, blood compatible fluid or CO₂.
2. If using CO₂, corrections to inflation volume may be necessary due to CO₂ diffusion through the balloon wall.
3. Purge the catheter and through-lumen of all air with a sterile medium.
4. For optimal control, the recommended syringe is included. See Specification Table.
5. Insert the catheter into the vessel either through an arteriotomy or appropriate introducer sheath, and position the balloon just beyond the occlusion.
6. Inflate the balloon with sterile medium. Inflation should continue until the balloon can be felt to engage the vessel wall. Sensitivity in controlling the balloon inflation volume during catheter withdrawal is important so that the balloon diameter can be adjusted. Inflation of the balloon is associated with a feeling of resistance against the vessel wall. When no resistance is encountered it should be assumed that the balloon has ruptured. Inflation should be discontinued immediately and the catheter withdrawn.
7. Withdraw the catheter to remove the occlusive material.

Complications

As with all catheterization procedures, complications may occur. These may include but are not limited to:

- Infection
- Local hematomas
- Intimal disruption
- Arterial dissection
- Perforation and rupture
- Hemorrhage
- Arterial thrombosis
- Distal emboli of blood clots or arteriosclerotic plaque
- Air embolus
- Aneurysms
- Arterial spasms
- Arteriovenous fistula formation
- Balloon rupture or tip separation with fragmentation and distal embolization

How Supplied

Syntel Over the Wire Embolectomy Catheters are supplied sterile in sealed peel-open packages. Sterility is assured as long as packages are unopened and/or undamaged. This device is supplied non-pyrogenic. Includes syringe.

Storage/Shelf Life

Handle with care. Product should be stored in a clean, cool, dry area away from chemical fumes.

Re-sterilization/Repackaging

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

Safe Handling and Disposal

This device is a single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

Cleaning

1. Devices considered necessary to return should be cleaned using one of the following:

- Sodium hypochlorite solution (500-600 mg/l), or
- Peracetic acid solution with subsequent ultrasonic treatment

2. Devices should then be decontaminated with either:

- 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
- Ethylene oxide gas

3. Devices should be completely dried prior to packaging.

Packaging

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.

2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.

3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.

4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.

5. Shipping paper and content marking for the outer shipping container are not required.

6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular
Attn: Complaint Lab
63 Second Avenue
Burlington, MA 01803, USA

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

| Model Number | Body Size | Maximum Liquid Capacity | Guidewire Compatibility | Deflated Compatible Introducer Sheath | Inflated Balloon Diameter (mm) | CO ₂ Volume (cc) | Max. Pull Force on Inflated Balloon (lbs.) | Recommended Syringe |
|--------------|-----------|-------------------------|-------------------------|---------------------------------------|--------------------------------|-----------------------------|--|---------------------|
| A4E01 | 5Fx40 cm | 0.15 ml | 0.035" | 6F | 5 mm ± 2.0 | 0.3 | 0.7 | 1cc |
| A4E02 | 5Fx80 cm | 0.15 ml | 0.035" | 6F | 5 mm ± 2.0 | 0.3 | 0.7 | 1cc |
| A4E03 | 5Fx40 cm | 0.5 ml | 0.035" | 6F | 9 mm ± 1.0 | 1.0 | 1.5 | 1cc |
| A4E04 | 5Fx80 cm | 0.5 ml | 0.035" | 6F | 9 mm ± 1.0 | 1.0 | 1.5 | 1cc |
| A4E05 | 5Fx40 cm | 0.9 ml | 0.035" | 6F | 11 mm ± 1.0 | 1.8 | 2.0 | 1cc |
| A4E06 | 5Fx80 cm | 0.9 ml | 0.035" | 6F | 11 mm ± 1.0 | 1.8 | 2.0 | 1cc |
| A4E08 | 5Fx80 cm | 1.4 ml | 0.035" | 6F | 13 mm ± 1.0 | 2.8 | 2.5 | 3cc |
| A4E09 | 5Fx80 cm | 1.6 ml | 0.035" | 6F | 14 mm ± 1.0 | 3.2 | 3.5 | 3cc |

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LeMaitre Vascular GmbH
LeMaitre Vascular Ltd