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Processed bovine pericardial patch

Instructions for use - English

STERILE **A**  **Rx only**

How to Use

A. Precautions

1. Carefully inspect the entire container and evidence of opening to check for any damage.
2. Do not use the product if the container is damaged or not sealed.
3. Clean surgical gloves sufficiently to remove all powder residues before handling this product.

B. How to Use

1. Selection of Product Model
 - 1) Select the appropriate type of model for the procedure which will be performed.
 - 2) Check the information on the container label to make sure that the correct size of the product is selected.
2. Cleaning procedure
 - 1) In order to reduce the patient exposure to residual glutaraldehyde, follow an appropriate cleaning procedure in accordance with the attached table. Clean several patches separately with new sterile saline solution.
 - 2) Remove the outer plastic seal and open the lid of the container. Make sure that the contents are sterilized, and treated in an aseptic condition to prevent contamination. The outer surface of the container is not sterilized, so do not move it to a sterile location.
 - 3) Remove the container while holding the edges of the product with the sterilized atraumatic forceps.
 - 4) Remove the product from the container, and soak it in sterilized saline solution. Gently stir the product in a cleaning basin using the same forceps. Leave the product in a cleaning basin until the surgeon asks for it.
 - 5) Because XenoSure bovine pericardial patch have been proven to have no side effects in the antibiotic treatment, depending on the permissions of the surgeon, the 500 U/mL of bacitracin or 10 mg/mL of cephalexin may be contained in a rinse solution. The effects of other antibiotics or the long-term effects of XenoSure bovine pericardial patch materials of these antibiotics have not been tested. Use antibiotics as directed by the manufacturer.
3. Cleaning procedure for alternate patch sizes
For the patch sizes not listed in the table below, follow the following cleaning instructions.
: Patch Length (cm) x Patch Width (cm) = Patch Area (cm²)
If the patch is smaller than 37.5 cm², clean it with 500 ml of saline solution for 2 minutes.
If the patch is larger than 37.5 cm² and smaller than 300 cm², clean it with 1000 ml of saline solution for 3 minutes. (The patch area is based on one side of the patch. This calculation is deliberately formalized. Follow this instruction.)

Table 1. Cleaning procedure according to each model

Model	Size (cm)	Cleaning procedure
0.6P8	0.6x8	500ml Clean in saline solution for 2 minutes
0.8P8	0.8x8	
1P6	1x6	
1P10	1x10	
1P14	1x14	
1.5P10	1.5x10	
2P9	2x9	

Model	Size (cm)	Cleaning procedure
4P4	4x4	500ml Clean in saline solution for 2 minutes
4P6	4x6	
6P8	6x8	1000ml Clean in saline solution for 3 minutes
8P14	8x14	
10P16	10x16	

4. Implantation

- 1) Cut or amend the product into the desired shape.
- 2) To prevent dryness during implantation, put the product in the sterilized normal saline solution frequently.
- 3) Visually inspect both sides of XenoSure® Biologic Patch. If one side appears smoother, implant the soft surface to contact bloodstream.
- 4) All surgeons who perform the above procedure must be properly trained and fully familiar with the relevant scientific paper.

C. How to Store after Use

1. This product is disposable so that it can not be reused or re-sterilized.
2. The leftover products must be classified as medical wastes and discarded accordingly.

Precautions for Use

A. Warning

The major complications reported for bovine pericardial tissue are fibrosis and infection. This complication is only observed for a few patients after the bovine pericardial implantation.

B. Caution

All medical personnel responsible for handling and preparing this product must pay particular attention to prevent the product damage.

- 1) This product is disposable so that it can not be reused or re-sterilized. Due to the failure of reuse, and/or re-process, and/or re-sterilization, it may result in patient injury, disease, or death. The unused parts of the product must be discarded. Check the expiration date of the product.
- 2) Check the sealed and sterilized containers before opening. Do not use the product if the packaging is damaged. Patients can be infected when the contents are not sterilized.
- 3) Clean the product according to the "Cleaning procedure" in the user manual before using the product. The product stock solution contains glutaraldehyde, so it can cause inflammation in skin, eyes, nose and throat. Do not inhale the vapour of the stock solution. Avoid prolonged skin contact and immediately wash the area with water. Medical measures must be taken immediately in case of contacting with eyes. Liquid chemical stock solutions must be discard in accordance with hospital procedures.
- 4) If the product is used with traumatic devices, it may damage the product.
- 5) Do not use the damaged products as the product integrity may have impaired.
- 6) If the product is damaged before implantation, replace it with another product without repairing it.
- 7) Do not re-sterilize the products. The unused parts are considered as non-sterilization so it must be discarded.
- 8) Do not expose the product to vapors, ethylene oxide, chemicals or radiation (gamma/electronic beam) sterilization that results in the product damage.
- 9) Do not use the cutting suture needles or the armed skin suture of the cutting point as it may damage the product.
- 10) Do not dry the product while handling it.
- 11) Do not use the expired products.

C. Side effects

This product is designed to restore the organs of the body. The inappropriate functions of the implanted product causes the same symptoms of a defect in organs. The implantation surgeon must inform the patients of any symptoms regarding the inappropriate functions of the product.

- 1) Complete heart block and right bundle branch blockage are the

reported complications for the procedures involving the heart recovery near A-V conductive bundle.

- 2) The bovine pericardium fixed with glutaraldehyde, might be attacked later through the immune system and can lead to the tissue deterioration. The benefits of using this product must be considered along with the risk possibility of late tissue damage.
- 3) Residual glutaraldehyde is at risk of toxic effects. When the appropriate rinse procedure mentioned in the user manual is completed, the residual glutaraldehyde in the patch will be reduced to the acceptable levels. Thus, the risk of acute toxic effects will be greatly reduced. In the review of the published paper, the safety limits for glutaraldehyde exposure are not established when it is implanted into bloodstream. In the case of implanting the tissues treated with a large amount of glutaraldehyde (e.g., multiple large patches), or for the patients with low weight, the risk will increase. The benefits of using this product must be considered along with the risk possibility of toxic effects.
- 4) Epithelial signs of calcification and malignancy were reported as side effects in the animal studies related to bovine pericardium. Due to the phagocytosis accompanied by the chronic inflammatory infiltration between the bovine pericardium and the surrounding host tissue, the implantation collagen degradation that is matched with the graft-versus-host reaction occurs.

D. Potential complications

- Repetitive stenosis
- Pseudoaneurysm formation
- Infection
- Thrombosis
- Calcification
- Fibrosis
- Vascular closure
- Patch rupture
- Expansion
- Myocardial infarction
- Bleeding
- Cerebrospinal fluid leak
- Stroke
- Death

Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and XenoSure are registered trademarks of LeMaitre Vascular, Inc.
© 2020 LeMaitre Vascular, Inc. All rights reserved.

EC REP

LeMaitre Vascular ULC
9135 Keele Street, Suite B6
Vaughan, Ontario
Canada L4K 0J4
Tel: 855-673-2266

LeMaitre Vascular GmbH
Otto-Volger-Str. 5a/b
65843 Sulzbach/Ts., Germany
Tel: +49-(0)6196-659230

LeMaitre Vascular, Inc.
63 Second Avenue
Burlington, MA 01803