**LeMaitre® Valvulotome**

Instructions for Use - English

**LeMaitre® Valvulotom**

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**Valvulotome LeMaitre®**

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**Valvulotomo LeMaitre®**

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**Valvulotomo LeMaitre®**

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**Valvulótomo LeMaitre®**

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**STERILE EO**

## Product Description

The LeMaitre® Valvulotome is a device that cuts venous valves during vascular procedures such as in-situ peripheral bypass and non-reversed translocated bypass. The Centering Hoops keep the head of the LeMaitre Valvulotome centered in the vessel and prevents the valve-cutting blades from damaging the vessel wall. The size of the Centering Hoops and Cutting Blades adjusts to the internal diameter of the vein as the LeMaitre Valvulotome is being drawn through the vessel.

The device contains a hydrophilic coating on the protective sheath and silicone coating on the metal components to improve trackability and reduce friction.

The LeMills Valvulotome is a device that cuts venous valves during vascular procedures. The device consists of a stainless steel blade and a plastic handle.

## Indication for Use

The LeMaitre Valvulotome and LeMills Valvulotome are indicated for the treatment of peripheral arterial disease which includes popliteal artery, femoral artery, tibioperoneal trunk, posterior tibial artery, peroneal artery, and peripheral vascular disease which includes atherosclerosis.

## Intended Use/Intended Purpose

The LeMaitre Valvulotome and LeMills Valvulotome are intended for excising or disrupting venous valves.

## Intended User

The LeMaitre Valvulotome and LeMills Valvulotome are surgical tools intended for use by experienced vascular surgeons trained in the procedures for which they are intended.

## Patient Population

Patients of any gender, age or ethnicity with vascular disorders. There is no data for the use of this device on pregnant women and children. It is at the surgeon's discretion on whether to use it on this population.

## Part of the Body Contacted

The LeMaitre Valvulotome and LeMills Valvulotome will come in contact with veins.

## Clinical Condition

In situ peripheral bypass, non-reversed translocated bypass.

## Clinical Benefits

The clinical benefits associated with the use of the LeMaitre Valvulotome and LeMills Valvulotome include increased limb salvage, amputation-free survival, and overall survival, comparable to rates observed with similar valvulotomes.

## Contraindications

1. Endarterectomy procedures
2. Thrombolysis procedures
3. Vein stripping procedures
4. Embolectomy procedures
5. Vessel dilation procedures

## Warnings

1. Do not use if the package or device is damaged.
2. Do not use if inner packaging is opened outside of a sterile environment.
3. The LeMaitre Valvulotome and LeMills Valvulotome are to be used by a qualified physician.
4. Do not insert the LeMaitre Valvulotome into a vessel or extract from a vessel in the open position.
5. Do not open or close the LeMaitre Valvulotome while in a coiled configuration.
6. The LeMaitre Valvulotome must be flushed with saline or heparinized saline only.
7. Do not pass the LeMaitre Valvulotome or LeMills Valvulotome through a vessel that has undergone synthetic grafting or contains implants.
8. Do not rotate the LeMaitre Valvulotome in a vessel.
9. Do not use device for valvulotomy unless the target vein is fully distended by arterial blood flow (due to anastomosis of artery to vein) or by saline injection.
10. Activate the hydrophilic coating with sterile saline before use. Do not use other solvents to activate the coating. Avoid wiping the device with dry gauze or other abrasive materials as this may damage the coating. Avoid excessive wiping of the coated device.
11. Do not advance or withdraw this device through sharp or beveled edges which may result in damage of the coating.
12. Failure to abide by these warnings may result in damage to the device coating which may necessitate intervention or result in a serious adverse event.

## Precautions

1. United States Federal and other law restricts this device to sale on or by the order of a physician.
2. Do not use these devices past the expiration date printed on the labeling.
3. These devices are single use only. Do not reuse, reprocess, or re-sterilize (reference "Re-sterilization/Repackaging" section of this document for further details).
4. Keep any dissected portion of the vein straight and avoid twisting the vein during valvulotomy.
5. Use caution when tracking the LeMaitre Valvulotome or the LeMills Valvulotome in thrombophlebitic veins.
6. Do not advance the LeMaitre Valvulotome with the blades in the open position.
7. Do not rotate the Green Handle while the LeMaitre Valvulotome device is in open position.

## List of Risks

### Potential device-related complications:


- Residual tributaries
- Entrapment of cutting heads/ blades in branches or anastomoses
- Inadequate valve disruption
- Vein stricture

**Potential procedure-related complications:**


- Vessel wall perforation
- Intimal disruption
- Thrombus formation
- Post-procedure wound necrosis
- Vessel damage
- Vessel occlusion/ stenosis
- Hematoma
- Hemorrhage
- Infections
- Erythema
- Embolism
- Death
- Hemolysis
- Irritation
- Toxicity

**Specifications**

LeMaitre Valvulotome:

LeMills Valvulotome:

Catalog Number	Useable Length	Closed Catheter Outer Diameter	Maximum Open Cutting Blade Diameter	Maximum Open Centering Hoop Diameter	LeMills Valvulotome	3cc Syringe
1009-00M	98 cm	1.5 mm	6.0 mm	9.5 mm	Included	Included
1010-00M	40 cm	1.5 mm	6.0 mm	9.5 mm	Not Included	Included

**To Open Package**

1. Open box and remove sealed sterile unit.
  2. Open sealed sterile unit and present contents to personnel in sterile environment.
  3. Carefully remove lids, as applicable.
- Caution:** Remove the LeMaitre Valvulotome from the tray with care to ensure centering hoops are not damaged.
4. Place in a basin filled with sterile, heparinized saline.

**Pre-Use Check**

Inspect blades for damage and alignment prior to use.

**Instructions for Use - LeMaitre Valvulotome**


1. Cut the first valve under direct vision using standard technique.
2. Perform the proximal anastomosis to distend the vein.
3. With the LeMaitre Valvulotome in the open position, remove the Irrigation Port cap and flush with sterile heparinized saline. Replace the cap upon completion.
4. Sheath the Centering Hoops within the Protective Sheath by advancing the White Handle until the Cutting Blades are fully enclosed and a positive stop is detected.
5. Activate the hydrophilic coating of the device by:
  - a) placing the LeMaitre Valvulotome into a basin with enough sterile saline to cover the device, or
  - b) wiping the Protective Sheath with sterile saline.
6. Advance the device into the vein to the treatment area; 2 cm to 3 cm distal to the anastomosis.
7. Expose the Cutting Blades and Centering Hoops by maintaining the position of the Green Handle and retracting the White Handle until a positive stop is detected.

**Note:** Failure to maintain the position of the Green Handle while exposing the Cutting Blades may cause damage to the anastomosis.

8. Slowly retract the device in the open position to disrupt the valve sets.

**Warning:** Do not rotate the LeMaitre Valvulotome.

9. Proceed with caution until the Green Safety Stripes appear. These stripes indicate the 2 cm end of the Protective Sheath. When the third (and final) Green Safety Stripe is visible, sheath the Cutting Blades and Centering Hoops by maintaining the position of the Green Handle and advancing the White Handle.

**Warning:** Failure to sheath the Cutting Blades and Centering Hoops prior to retracting the device from the vein may cause damage.

**Note:** If sheathing the Cutting Blades for vein extraction leaves the most distal valve set intact, use the LeMills Valvulotome to excise using standard technique.

**Note:** If the device is set aside for further use within the same procedure, it must be stored in the open position in a basin of heparinized saline. Flush the device as noted in step 3 prior to reintroduction.

10. Pass the device through the vein an additional 1-2 times.
11. Steps 5 through 9 may be repeated as necessary.
12. Confirm the free flow of blood through the vein.

### Instructions for Use - LeMills Valvulotome

1. Insert the device into the vein to the treatment area.
2. Slowly retract the device to disrupt the valve sets.
3. Steps 1 and 2 may be repeated as necessary.
4. Confirm the free flow of blood through the vein.

### Storage

Store in a cool, dry place.

### Resterilization/Re-use

These devices are single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

### Safe Handling and Disposal

This device is single-use and disposable device. Do not implant. Please return the used device only at the time that the device has not performed as intended or the device is related to an adverse event. In other situations, the device should not be returned but disposed according to local regulations.

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no heavy metals or radioisotopes, and is not infectious or pathogenic. Please consult local regulations to verify proper disposal of sharps.

### Cleaning:

1. Devices considered necessary to return should be cleaned using one of the following:
  - a) Sodium hypochlorite solution (500-600 mg/l), or
  - b) Peracetic acid solution with subsequent ultrasonic treatment
2. Devices should then be decontaminated with either:
  - a) 70% solutions of ethanol or isopropanol for a minimum of 3 hours or,
  - b) Ethylene oxide gas
3. Devices should be completely dried prior to packaging.

### Packaging:

1. Cleaned devices should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. For devices capable of penetrating or cutting skin or packaging material, the primary packaging must be capable of maintaining the product without puncture of the packaging under normal conditions of transport.
2. The sealed primary container should be placed inside watertight secondary packaging. The secondary packaging should be labelled with an itemized list of the contents of the primary receptacle. Cleaning methods should be detailed if possible.
3. Both primary and secondary packaging of cleaned, decontaminated single-use disposable devices should be labelled with an ISO 7000-0659 Biohazard symbol.
4. Primary and secondary packaging must then be packaged inside an outer package, which must be a rigid, fiberboard box. The outer shipping container must be provided with sufficient cushioning material to prevent movement between the secondary and outer containers.
5. Shipping paper and content marking for the outer shipping container are not required.
6. Packages prepared in the above manner may be shipped to:

LeMaitre Vascular  
Attn: Complaint Lab  
63 Second Avenue  
Burlington, MA 01803, USA

### Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc., warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME. This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre Vascular GmbH
LeMaitre Vascular Ltd