CardioCel™ Bioscaffold Patch

Patient Information – English

CardioCel™ Biogerüstpatch

Patienteninformation – Deutsch

Patch Bioscaffold CardioCel™

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CardioCel Neo Bioscaffold Patch	Width	Length	Model #
Collagen Bioscaffold Patch	4 cm	4 cm	EC0404N
Collagen Bioscaffold Patch	5 cm	8 cm	EC0508N

Indication for Use

The cardiac patch is indicated for use in the repair of heart and blood vessels.

Intended Use

The cardiac patch is used in heart and blood vessel defects. The patch material is an implant to repair damaged heart or blood vessel tissue.

Intended User

The users of the patch are qualified surgeons.

Patient Population

Patients of any gender, age or ethnicity in need of an implant to repair heart deformities. There is no data for the use of this device on pregnant women.

Part of the Body Contacted

- Heart tissue and blood vessels

Clinical Conditions

- Cardiac and vascular defects

Clinical Benefits

- Increased survival rates
- Improved quality of life:
 - General improvement in overall health/ wellness
 - Improvement in exercise tolerance
- Prevention/reduction of further surgery in later life

Self-care Instructions

1. Your new device is a foreign body and therefore needs close monitoring and careful observation. It may take 6-8 weeks for full recovery.
2. After placement, the implant area may be swollen and tender for up to a week.
3. Watch for any new redness or tenderness.
4. Watch for any opening in the incision(s).
5. Watch for numbness, tingling or pain.

NOTE: If you experience any symptoms described in 3, 4 or 5 above please contact your provider.

6. Do not puncture or manipulate the patch.
7. If the patch was implanted in your leg, swelling in the extremity is expected because of increased blood flow. Elevate or move the extremity according to your provider's instructions.
8. It is preferable to have the surgical site covered for the first week to protect skin and incision(s). (Follow your provider's instructions.)
9. Keep bandages or wound covering on as per your provider's instructions.
10. If you have adhesive surgical tape or strips across your incision(s), wear loose clothing that does not rub against your incision(s). The adhesive surgical tape or strips will curl up and fall off on their own after a week.
11. You may shower or get the incision(s) wet, once your provider says you can. DO NOT soak, scrub, or have the shower beat directly on the incision(s).
12. DO NOT soak in the bathtub, a hot tub, or a swimming pool. Ask your provider when you can start doing these activities again.
13. Your provider will tell you how often to change your wound covering and when you may stop using one. Keep your incision(s) dry. If your incision(s) goes to your groin, keep a dry gauze pad over it to keep it dry.
14. Clean your incision(s) with soap and water every day once your provider says you can. Look carefully for any changes. Gently pat it dry.
15. DO NOT put any lotion, cream, or herbal remedy on your incision(s) without first discussing with your provider.
16. Consult your provider for instructions on taking any prescription or over-the-counter medications after surgery.

Intended Performance

- The patch is intended to repair heart and blood vessel defects.

Undesirable Side Effects

Potential device-related adverse events

- When an artery that was opened becomes narrowed again (Restenosis)
- A life-threatening inflammation of the inner lining of the heart's chambers and valves (Infective endocarditis)
- Buildup of excess calcium (Calcification)
- Ruptured red blood cells (Haemolysis)
- Blood clots in veins (Thromboembolism)
- Inflammation
- Decline (degeneration) of the implants
- Formation of clinically significant fibrous tissue
- Infection
- Blood clotting in vein (Thrombosis)
- Graft becomes dilated (Dilatation)
- Heart Attack (Myocardial infarction)
- Bleeding

- Stroke
 - Death
- Potential procedure-related adverse events*
- Narrowing of tubular structures (Stenosis)
 - Flow obstruction
 - When your heart can't pump well because the pouch around it thickens (Pericardial adhesions)
 - Inflammation
 - Injured blood vessel wall that leads to leaking (Pseudoaneurysm formation)
 - Patch rupture

Long-term Protection Measures for Your Patch

- For cardiac implants there are no long term protection measures. Please contact your provider for further guidance.
- If implanted in an extremity:
 - avoid prolonged extreme extension of the arm or leg with the implantation as it could lead to nerve damage.
 - avoid extreme or abrupt movements of the arm, shoulder, or legs during a post-operative period of 1.5 to 2 months. Specifically, you should not reach out in front, raise arms above shoulder level, throw, pull, stride, or twist.
 - avoid sleeping on the patch implantation side of your body as it may cause compression.

Risks Related to Interactions with Other Equipment

- None

Post-surgical Monitoring

- Check your incision(s) every day.
- Call your provider immediately if you have any signs of a blood clot, swelling, unusual skin color or infection, such as:
 1. Increased pain
 2. Swelling, redness, or red streaks
 3. Blood or pus draining from the incision(s)
 4. Numbness
 5. Fever
- The nature and frequency of regular or preventative examination, monitoring or maintenance of the patch will be determined by your provider. This will be based on your underlying medical condition and the status of your patch.

Lifetime of the Patch

The lifetime data of the patch is based on the conditions below:

- Repair of holes/ abnormalities in the heart, up to 10 years
- Procedure to treat heart valves and heart tissue, up to 5 years
- Repair of major blood vessels that connect to the heart, up to 10 years

Proper care and regular medical follow up may extend the life of the device.

To ensure your patch functions as intended, follow the guidance of your healthcare provider.

When to Contact Your Provider

- Closely watch for any changes in your health. Seek emergency care anytime you experience:
 - Lost consciousness
 - Trouble breathing
 - Extremity has severe pain or becomes cold, pale, blue, tingly, or numb
 - Pain that does not get better after you take pain medicine
 - Loose stitches or your incision(s) opens
 - Extensive bleeding from the incision(s)
 - Signs of infection, such as:
 1. Increased pain, swelling, warmth, or redness
 2. Red streaks leading from the incision(s)
 3. Pus draining from the incision(s)
 4. A fever over 101 °F (38.3°C)
 5. You are nauseous or cannot keep fluids down
 - You have chest pain, dizziness, problems thinking clearly, or shortness of breath that does not go away when you rest
 - You are coughing up blood or yellow or green mucus
 - You have chills
 - You experience abdominal pain or bloating

Device Materials

The patient may be exposed to the following materials and substances for worst-case largest patch 6 cm x 14 cm:

- Processed tissue from cow hearts: average 5.55g per patch
- Propylene glycol: up to 13.0mg per patch
- Propylene oxide: less than detectable limit
- Glutaraldehyde: less than detectable limit

The CardioCel patch has passed testing to ensure its biocompatibility.

In Case of an Emergency

- Any serious adverse reaction should immediately be reported to your provider.
- Any device related emergency that occurs in the EU should be reported to the Competent Authority of the country where you are located and the manufacturer, LeMaitre Vascular, Inc. (+1 781-221-2266 or; <https://www.lemaitre.com/contact-us/email-us>).
- Any device related emergency that occurs in Australia should be reported to the Therapeutic Goods Administration via their website (<https://www.tga.gov.au/>) and the manufacturer, LeMaitre Vascular, Inc. (+1 781-221-2266 or; <https://www.lemaitre.com/contact-us/email-us>).

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LeMaitre Vascular GmbH
LeMaitre Vascular, Ltd.