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AlboSure® Polyester Vascular Patch

AlboSure® Polyester Vascular Patch

Instructions for Use – English

STERILE EO

Description

AlboSure® Vascular Patches, made of a synthetic polyester (polyethylene terephthalate, PET) material, are designed to repair and/or close the cardiovascular system. The fabric is designed with a run-proof structure reducing the risk of fraying or wearing down at edges.

A guideline facilitates implantation.

The AlboSure is impregnated with collagen, reducing permeability and mitigating hemorrhaging problems during implantation. This eliminates the need to preclot prior to the procedure. The process of impregnation using bovine collagen maintains the original structure of the material and the structural characteristics of the prosthesis, i.e. flexibility and softness. The collagen is cross-linked to the graft by formaldehyde steaming, ensuring a balanced reabsorption rate by the host tissues.

Indication for Use

The AlboSure Polyester Vascular Patch is indicated for cardiac and vascular patch grafting. The fabric is recommended for use in patients requiring systemic heparinization prior to, or during surgery.

Permeability

Permeability is a chief property of the polyester fabric and depends on the porosity of the material. This capability is expressed as the amount of water that can pass through 1 cm² of the prosthesis at 120 mmHg over a given unit of time. The permeability of AlboSure is $\leq 1 \text{ ml cm}^2 \text{ min}^{-1}$.

Packaging

The AlboSure Vascular Patch is supplied sterile. The product is supplied in a double pouched package; however, the inner pouch is internally sterile only. The inner product sleeve is externally sterile, and can be brought into sterile field. The outer pouch is not sterile; it must not be brought into a sterile environment.

Sterility

The AlboSure Vascular Patch contained in this package ("this device") is supplied sterile unless the inner package is opened or damaged. Be careful to inspect this device and its packaging to verify that no damage has occurred as a result of shipping, and that the sterile barrier has not been compromised. The AlboSure Vascular Patch must never be resterilized.

Instruction For Use

Choice and Preparation

Choose the patch size best suited to the vascular section to be replaced. Trimming the patch to proper size is critical, and should be done using sharp surgical scissors. Carefully examine the packaging enclosing the selected patch to verify its integrity, the absence of holes, cuts or perforations, and the expiration date.

Remove the patch from package using an aseptic procedure.

The collagen impregnated vascular patches can be used as supplied, or after immersing in a sterile saline solution just before implantation, making it flexible and easy to handle. During the preparation and procedure, care should be taken not to put excessive tension on the patch. Correct handling helps avoid the risk of damaging the collagen impregnation on the patch surface.

Surgical Technique

Considering how complex a surgical implantation of a vascular patch is, and how many factors can influence the patient's condition, it is left to the surgeon to define suture technique and patch size, as well as the therapy to adopt before, during and after the procedure.

Contraindications

The patch is not indicated for use in:

- Reconstruction of hernias or tissue deficiencies.
- Reconstruction or repair of passive biological membranes (such as dura mater, pericardium, or peritoneum).

Precautions and Warnings

- Do not preclot.
- Do not attempt to re-sterilize or reuse this device. Do not use if the package is damaged or opened outside a sterile field.
- Avoid contact with extraneous particles (contact may generate emboli or undesirable interactions with the blood).
- Surgical gloves used to handle patch should not contain powders, preservatives or lubricants.
- Avoid damaging the patch when handling. Atraumatic clamps and appropriate instruments (e.g. vascular clamps) should be used when handling is required. Avoid using these instruments with undue force, otherwise the collagen coating or fabric will be damaged.
- Avoid excessive tension on fabric.
- Atraumatic needles are recommended.
- Sharp surgical scissors are recommended to trim patch to proper size.
- The patch should not be implanted in patients who exhibit sensitivity to polyester or materials of bovine origin.

Adverse Effects

- | | |
|----------------------------|-----------------------------|
| • Restenosis | • Dilatation |
| • Arterial rupture | • Myocardial infarction |
| • Pseudoaneurysm formation | • Bleeding |
| • Infection | • Stroke |
| • Thrombosis | • Death |
| • Vessel occlusion | • Cranial Nerve Injury |
| • Patch rupture | • Transient Ischemic Attack |

Storage

Store in a cool, dry place at a temperature between 0°C and 30°C.

Resterilization/Re-use

This device is single-use only. Do not reuse, reprocess, or re-sterilize. The cleanliness and sterility of the re-processed device cannot be assured. Reuse of the device may lead to cross contamination, infection, or patient death. The performance characteristics of the device may be compromised due to reprocessing or re-sterilization since the device was only designed and tested for single use. The shelf life of the device is based on single use only. If for any reason this device must be returned to LeMaitre Vascular, place it in its original packaging and return it to the address listed on the box.

Safe Handling and Disposal

If serious medical incidents should arise during use of this medical device, users should notify both LeMaitre Vascular and the Competent Authority of the country where the user is located.

This product contains no sharps, heavy metals or radioisotopes, and is not infectious or pathogenic. No special requirements for disposal are evident. Please consult local regulations to verify proper disposal.

Packaging and Shipping of Explanted AlboSure:

Return of the shipment to LeMaitre Vascular depends on 3 crucial questions:

1. Is the explant obtained from a patient with a known or presumed pathogenic condition at the time of the explant?
2. Is the explant obtained from a patient that has a known treatment history that involves therapeutic radionuclides within the last 6 months?
3. Has the clinician obtained consent from the patient for the sample to be returned to the manufacturer for research purposes?

In the event that the response to question 1 or 2 is affirmative, LeMaitre Vascular does not provide adequate guidance for shipment. THESE EXPLANTS SHOULD NOT BE RETURNED TO LEMAITRE VASCULAR UNDER ANY CIRCUMSTANCES. In these cases, the explant should be disposed according to local regulations.

For explants that do not have pathogenic or radiologic hazards, please use the following:

Pre-explantation:

1. If possible, perform a CT or Ultrasound scan of the device to document patency.
2. LeMaitre Vascular can accept clinical information that is patient-anonymized. LeMaitre Vascular requests

information including:

- a) The original diagnosis which resulted in use of the implant,
- b) The patient's medical history relevant to the implant, including the hospital or clinic at which the device was implanted.
- c) The patient's implant experience prior to implant removal.
- d) The hospital or clinic at which the explantation was done, and date of retrieval.

Explantation:

1. Explanted AlboSure patches should be directly transferred to a sealable container filled with a solution of alkaline buffered 2% glutaraldehyde or 4% formaldehyde prior to shipment.
2. Cleaning of explanted patches should be minimal if necessary. Proteolytic digestion should not be used in any circumstances.
3. AlboSure explants should not be decontaminated under any circumstances. DO NOT autoclave the sample or use ethylene oxide gas to decontaminate.

Packaging:

1. Explants should be sealed and packed in a manner that minimizes potential for breakage, contamination of the environment or exposure to those handling such packages during transit. Material that is absorbent and is cushioning should be selected for isolating the sealable container inside the secondary packaging. Primary and secondary packaging must then be packaged inside an outer package.
2. Explants in sealed primary containers should be labeled with an ISO 7000-0659 Biohazard symbol. The same symbol should be attached to the secondary packaging and to the outer packaging. Outer packaging should also be labelled with Name, Address and Telephone Number of Sender, and the statement, "Upon discovery of damage or leakage, the package should be isolated and sender notified".
3. Packages prepared in the above manner may be shipped to:
LeMaitre Vascular
63 Second Avenue
Burlington, MA 01803

Notices: Limited Product Warranty; Limitation of Remedies

LeMaitre Vascular, Inc. warrants that reasonable care has been used in the manufacture of this device and that this device is suitable for the indication(s) expressly specified in these instructions for use. Except as explicitly provided herein, LEMAITRE VASCULAR (AS USED IN THIS SECTION, SUCH TERM INCLUDES LEMAITRE VASCULAR, INC., ITS AFFILIATES, AND THEIR RESPECTIVE EMPLOYEES, OFFICERS, DIRECTORS, MANAGERS, AND AGENTS) MAKES NO EXPRESS OR IMPLIED WARRANTIES WITH RESPECT TO THIS DEVICE, WHETHER ARISING BY OPERATION OF LAW OR OTHERWISE (INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND HEREBY DISCLAIMS THE SAME.

This limited warranty does not apply to the extent of any abuse or misuse of, or failure to properly store, this device by the purchaser or any third party. The sole remedy for a breach of this limited warranty shall be replacement of, or refund of the purchase price for, this device (at LeMaitre Vascular's sole option) following the purchaser's return of the device to LeMaitre Vascular. This warranty shall terminate on the expiration date for this device.

IN NO EVENT SHALL LEMAITRE VASCULAR BE LIABLE FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, SPECIAL, PUNITIVE, OR EXEMPLARY DAMAGES. IN NO EVENT WILL THE AGGREGATE LIABILITY OF LEMAITRE VASCULAR WITH RESPECT TO THIS DEVICE, HOWEVER ARISING, UNDER ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE, EXCEED ONE THOUSAND DOLLARS (US\$1,000), REGARDLESS OF WHETHER LEMAITRE VASCULAR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, AND NOTWITHSTANDING THE FAILURE OF THE ESSENTIAL PURPOSE OF ANY REMEDY. THESE LIMITATIONS APPLY TO ANY THIRD-PARTY CLAIMS.

These limitations do not apply to consumers in Australia or to the extent they are precluded by local law in any other jurisdiction.

A revision or issue date for these instructions is included on the back page of these Instructions for Use for the user's information. If twenty-four (24) months has elapsed between this date and product use, the user should contact LeMaitre Vascular to see if additional product information is available.

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LeMaitre and AlboSure are registered trademarks of LeMaitre Vascular, Inc.
©2020 LeMaitre Vascular, Inc. All rights reserved. Printed in the USA.
LeMaitre Vascular GmbH